MicroMega One Glider

Reference: 30006062-B

	∅	%	L	Stop	Ring	RPM	Maximum Recommended Torque (N.cm)
MicroMega One Glider	16	03	21 - 25 - 31	Yellow	White	250-400 rpm	1.2

Instruction for use EN

1. Indication
Glide path of root canals during non-surgical endodontic treatment. For use by a dental professional only.

2. Contraindications
Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there are no contraindications to the use of MicroMega One Glider for endodontically treating a tooth by orthograde route.

3. Complications
In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

- Warnings/precautions**
 - The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
 - Patients identified as having a risk of infectious endocarditis.
 - Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.
 - Respect the good dental practice in particular the use of a dental dam and gloves.
 - Use in continuous rotation at the recommended speed.
 - In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.
 - Use according to the recommended procedure (§7).
 - Instruments supplied sterile! Check the integrity of the packaging before use. If damaged, do not use the instrument.
 - If applicable, respect the expiry date.
 - Keep the UDI information which is only on the labelling until the last use.
 - Year of manufacture: see labelling.
 - Do not use if doubt on the respect of storage conditions.
 - Single use instrument: Do not reuse the instrument. If you do, there is a risk of decreased performance and security characteristics (fracture of the instrument in the tooth).
 - In case of doubt concerning the product identification, do not use.
 - Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
 - Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

Packaging symbol	STERILEEO		Sterilized using Ethylene Oxide. Do not re-use
------------------	------------------	---	--

- Clinical claims**
MicroMega One Glider in normal conditions.

Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

Usability:

- Secure progression and limit the interference

Safety:

- Limited risk of breakage
- Limited risk of shaping aberration
- Limited risk of screwing

- Characteristics**
Type 4 rotating instruments under EN ISO 3630-1:2019 standard:
 - Maximum recommended torque (not required): see table «Recommended Torque».

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (Type 1).
Number of uses: single use. 4 canals maximum recommended on the same patient
Medical device class according to Directive 93/42/EEC and MDR 2017/745: IIa.

- Protocol**
General information on the protocol for use:

After opening the access cavity, sodium hypochlorite must be constantly used and renewed before and after each instrument. The sodium hypochlorite must be placed using a syringe in order to be injected as close as possible to the apical region.

Operating procedure :
MicroMega One Glider is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.

- Reprocessing instructions**
Not applicable.

- Storage and transport conditions**

 - Do not use if the instrument is damaged or shows signs of wear, do not use it.

- Disposal**
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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