EASY CANALS		
Diameter	6%	4%
N.30	2/3 WL	WL

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CANALS OF AVERAGE DIFFICULTY		
Diameter	6%	4%
N.25	2/3 WL	WL
N.30		WL

DIFFICULT CANALS		
Diameter	6%	4%
N.20	2/3 WL	WL
N.25		WL
N.30		WL

Diameter	%	L	Stop	Ring	RPM	Recommended torque (N.cm)
20	4	21	Grey	Yellow	300-600	0.6 - 1
		25				
		29				
		6				
25	4	21	Grey	Red		0.6 -1
		25				
		29				
		6				
30	4	21	Grey	Blue		1 - 1.5
		25				
		29				
		6				

1. Indication
- Shaping of root canal during non-surgical endodontics treatment.
- For use by a dental professionals only.
2. Contraindications
- There is no contraindication to the use of HERO Shaper files for endodontically treating a tooth by orthograde way.
3. Complications
- In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.
4. Warnings/precautions
- The decision to use an endodontics instrument must be relayed to the clinical case expertise, particularly where the canal anatomy is considered to be too complex.

Respect recommendations for patients identified as having a risk of infectious endocarditis.

Contains Nickel and Titanium and must not be used on patients with known allergic sensitivity to these metals.

Respect the good dental practice in particular by using a dental dam and gloves.

Use in continuous rotation at the recommended speed.

In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.

Use according to the recommended procedure (S7).

Do not use in endodontics retreatment.

Respect the expiry date.

Keep the UDI information which is only on the labelling until the last use.

Year of manufacture : see labelling.

Do not use if you doubt on the respect of storage conditions.

In case of doubt concerning the product identification, do not use.

Check the condition of the instrument before the use between each canal. If the instrument is damaged or shows signs of wear, do not use it.

Inform the manufacturer and the national regulatory authority of any serious incident related to the instrument.

Be careful to respect the recommended maximum number of uses.

Instrument supplied non-sterile: follow the reprocessing procedure (S8) before use.
- | | | |
|------------------|--|-------------|
| Packaging symbol | | Non-sterile |
|------------------|--|-------------|
5. Clinical claims

Clinical performance in normal conditions of use :

- Shaping ability: adapted debris removal and cutting efficiency

- Respect of canal anatomy (respect of the original canal axis)

- Geometry adapted to irrigation

- Maintain of the apical foramen in its original position

Safety performance in normal conditions of use :

- Limited risk of breakage

- Limited risk of shaping aberration

- Limited risk of screwing

6. Characteristics

Type 2 continuous rotating instruments under EN ISO 3630-1: 2019 standard:

Maximum recommended torque: see table 'Recommended Torque'.

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (type 1).

Number of uses : maximum of 5 canals recommended if the instrument is not visually damaged.

Medical device class according to Directive 93/42 and MDR 2017/745: Ila

7. Protocol

EWL : Estimated Working Length

Initial penetration:

The first step consists of the initial penetration into the root canal with a conventional stainless steel hand instrument (usually a K file no.010) which provides information on the root canal anatomy complementing the information obtained through to the preoperative x-ray.

The mineralization of the root canal opening and the degree of curvature determine the difficulty of the treatment and thus the choice of the best suited sequence for the clinical case.

Penetration depths:

The .06 taper instruments are brought down to the 2/3 of the working length (WL). The .04 taper instruments are brought down to the WL. It is recommended to precisely determine the working length between the passage of the .06 taper instrument and the .04 taper instrument.

Operating dynamics:

The HERO SHAPER instrument is positioned at the canal entrance. Start the motor and maintain the speed at a constant level during the instrument's work until it is pulled out of the root canal.
- The instrument is brought to the desired penetration depth with an up and down movement of low amplitude with slight pressure applied to the contra-angle head. It is then possible to prepare the root canal walls with a circumferential brushing movement allowing to respect the root canal anatomy and prepare the lateral zones.

Irrigation:

It is important to abundantly irrigate the root canal with sodium hypochlorite.

Three sequences are available depending on the root canal difficulty.

Yellow sequence:

Sequence for difficult, strongly curved root canals possibly with intense mineralization of the root canal opening, making the initial penetration tricky, even with thin diameter limits.

Subsequent use of HERO Shaper® no.020 6%, no.020 4%, no.025 4% and no.030 4%.

Red sequence:

Sequence for root canals of medium difficulty which may present a stronger canal curvature; the passage of files for glide path development no. 10 down to the apex may be difficult.

Subsequent use of HERO Shaper® no.025 6%, no.025 4% and no.030 4%.

Blue sequence:

Sequence for root canals with a simple preparation, shallow curvature and a root canal opening which is sufficiently large for the passage of a file for glide path development no. 15 down to the apex in a first-line treatment.

Subsequent use of HERO Shaper® no.030 6% and no.030 4%.
- | 8. Reprocessing instructions | |
|------------------------------|--|
| General recommendations | <div><div>For your own safety, please wear personal protective equipment (gloves, glasses and mask).</div><div>Do not use cleaning agents containing phenol, aldehyde and alkaline composition.</div><div>Always respect the instruction for use provided by the product manufacturer.</div></div> |
| Limitations on preprocessing | <div><div>Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of medical devices is determined by their function and careful handling. Multiple use of cleaning and re-sterilization cycles may lead to increased risk of file separation.</div><div>The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements.</div><div>The state of the art and national laws require that validated procedures be followed.</div></div> |
| Material needed | <div><div>Gloves, masks, blouse, dental dam, as recommended by the cleaning agents manufacturers.</div><div>Running Tap water or deionized water.</div><div>Cleaning agent (neodisher® Septo Active and neodisher® MediZym)</div><div>Small soft brushes</div><div>Container</div><div>Ultrasonic tub or washer-disinfector</div><div>Class B sterilization apparatus according to EN ISO 17665</div></div> <div>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</div> |
- | | | |
|-----------|-----------------------------|--|
| 1 | Initial treatment | <div><div>Place used products in a container or a wipe with running tap water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min</div><div>Rinse the products with running tap water at 20-40 °C for 1 min.</div><div>If there is time before the next step, make sure the device remains moist by packing it in a wet wipe. Do not exceed 1 hour.</div></div> <div>Remarks:</div> <div>Do not use fixating agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning.</div> <div>Follow the instructions and respect the concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects to devices).</div> |
| | | |
| 2 | Preparation before cleaning | <div><div>If there are visible impurities on a device, it is recommended to manually pre-clean by brushing under running tap water at 20-40°C for at least 1 min with a soft brush until the impurities have been removed.</div></div> <div>Remark: Follow the instructions and respect the concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects to devices).</div> |
| | | |
| 3 | Visual inspection | <div><div>Inspect the used products and discard the damaged ones (broken, untwisted or anomalously bent).</div></div> |
| | | |
| manual | 4 | <div><div>Cleaning</div><div>Insert the products into an ultrasound apparatus beaker.</div><div>Run ultrasound apparatus during 10-15 min with deionized water and 0,5-2,0% neodisher® MediZym.</div></div> <div>Remarks: Follow the instructions, observe the water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution</div> |
| | | |
| | | <div><div>Rinsing</div><div>Rinse the products with tap water at 20-40°C for 1 min.</div></div> |
| | | |
| | | <div><div>Drying</div><div>Dry the products with compressed air until products are visibly dried.</div></div> |
| | | |
| automatic | 4 | <div><div>Cleaning/Rinsing/Drying</div><div>Place the instruments in the tray of the on the slide-in trolley of the washer/disinfector. Ensure that the cleaning agent and water can reach all area of the instrument's blade.</div><div>Pre-cleaning: cold tap water, 4 min</div><div>Run cleaning cycle with 0,5%-1,0% neodisher® MediZym at 50°C for 10 min.</div><div>Final rinsing with deionized water for 3 min.</div><div>Perform drying.</div></div> <div>Remarks:</div> <div>Follow the instructions and concentrations given by the detergent manufacturer.</div> <div>Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer.</div> <div>Use only approved washer-disinfector according to EN ISO 15883, maintained and validated regularly.</div> |
| | | |
| | | |
| | | |
| | | |
| | | |
| 7 | Visual inspection | <div><div>Inspect the used products.</div><div>Repeat steps 4, 5 and 6 if the product is not visibly clean or discard any damaged products.</div></div> |
| | | |
| 8 | Packaging | <div><div>Place the instruments twice in paper-plastic pouches for steam sterilization in compliance with ISO 11607 and EN 868.</div></div> <div>Remarks:</div> <div>Seal the pouches according to the manufacturer's recommendations. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.</div> |
| | | |
| 9 | Sterilisation | <div><div>Sterilize the products with steam vapor</div><div>Apparatus: class B</div><div>Temperature: 134°C</div><div>Time: 18min</div><div>Minimum drying: 20min</div><div>Check the physico-chemical indicators and cycle parameters.</div></div> <div>Remarks:</div> <div>When sterilizing multiple instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded.</div> <div>Place the pouches in the steam sterilizer according to the recommendation given by the sterilizer manufacturer.</div> <div>Use only pre-vacuum air removal steam sterilizer that are matching the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full size sterilizer), with saturated steam.</div> <div>Use only validated sterilizer according to EN ISO 17665, maintained and validated regularly</div> |
| | | |
| 10 | Storage | <div><div>Store the products in a dry, clean and dust-free environment at the temperature specified on the paper-plastic pouch by the steam sterilization manufacturer.</div></div> <div>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date). In case of damage, a complete rework should be performed.</div> |
| | | |
- | | | | |
|--|---|--|--------------------------|
| 9. Storage and transport conditions | | | |
| | | | |
| 10. Disposal | | | |
| After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.. | | | |
| 11. Symbols | | | |
| | NiTi raw material | | Root canal preparation |
| | Quantity | | Assortment |
| | EN ISO 1797 Type 1 | | Medical device |
| | Do not use if the packaging is damaged | | Legal manufacturer |
| | Sterilizable in a steam sterilizer (autoclave) at the specified temperature | | Unique Device Identifier |
- Year of CE-marking : 2001
Revision date : 2024-03-26
Reference : 80000231-A

0459

Micro-Mega SA
12, Rue du Tunnel
25000 Besançon
France
T +33 3 81 54 42 42
F +33 3 81 54 42 30
commercial@micro-mega.com
www.micro-mega.com

MicroMega