
One G NiTi Root Canal Instruments

Instruction for use EN

1. Indication
Glide path of root canals during non-surgical endodontic treatment. For use by a dental professional only.

2. Contraindications
Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there are no contraindications to the use of One G for endodontically treating a tooth by orthograde route.

3. Complications
In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc) could occur and lead to a risk of infectious processes.

- 4. Warnings/precautions**
- The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
 - Patients identified as having a risk of infectious endocarditis.
 - Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.
 - Respect the good dental practice in particular the use of a dental dam and gloves.
 - Use in continuous rotation at the recommended speed.
 - In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.
 - Use according to the procedure recommended (§7).
 - Instruments supplied sterile: Check the integrity of the packaging before use. If damaged, do not use the instrument.
 - If applicable, respect the expiry date.
 - Keep the UDI information which is only on the labelling until the last use.
 - Year of manufacture: see labelling.
 - Do not use if you doubt on the respect of storage conditions.
 - Single use instrument: Do not reuse the instrument. If you do, there is a risk of decreased performance and security characteristics (fracture of the instrument in the tooth).
 - In case of doubt concerning the product identification, do not use.
 - Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
 - Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

Packaging symbol			Sterilized using Ethylene Oxide. Do not re-use
------------------	--	--	--

5. Clinical claims
One G in normal conditions.

Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Maintains the apical foramen in its original position

Usability:

- Secure progression and limit the interference

Safety:

- Limited risk of breakage
- Limited screwing effects

6. Characteristics
Type 2 rotating instruments under EN ISO 3630-1: 2019 standard:

- Maximum recommended torque (not required): see table 'Recommended Torque'.

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (Type 1).

Number of uses: single use.

Medical device class according to Directive 93/42/EEC and MDR 2017/745: Ila.

7. Protocol Operating procedure :

One G is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.

8. Reprocessing instructions

Not applicable.

9. Storage and transport conditions

10. Disposal
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

11. Symbols

	NiTi raw material		Root canal preparation
	Quantity		Assortment
	'CLASSICS' handle		Do not use if packaging is damaged
	Medical device		

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