HERO642

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FR

EASY CANALS			
Diameter	6%	4%	2%
N. 30	1/2 - 2/3 WL	WL minus 2mm	WL

DIFFICULT CANALS			
Diameter	6%	4%	2%
N. 20	1/2 - 2/3 WL	WL minus 2mm	WL
N. 25		WL minus X	WL
N. 30			WL

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CANALS OF MEDIUM DIFFICULTY			
Diameter	6%	4%	2%
N. 25	1/2 - 2/3 WL	WL minus 2mm	WL
N. 30		WL minus X	WL

Diameter	%	L	Stop	Ring	RPM	Recommended torque (N.cm)
20	2	21	White	Yellow	300-600	0.6 - 1
		25				0.6 - 1
	4	21	Grey	1 - 1.5		
		25		0.6 - 1		
	6	21	Black	1 - 1.5		
		25		0.6 - 1		
25	2	21	White	Red	300-600	0.6 - 1
		25				0.6 - 1
	4	21	Grey	1 - 1.5		
		25		0.6 - 1		
	6	21	Black	1 - 1.5		
		25		0.6 - 1		
30	2	21	White	Blue	300-600	0.6 - 1
		25				1 - 1.5
	4	21	Grey	1 - 1.5		
		25		1.5 - 2		
	6	21	Black	1.5 - 2		
		25				

1. Indication
Shaping of root canal during non-surgical endodontics treatment.
For use by a dental professionals only.

2. Contraindications
There is no contraindication to the use of HERO 642 files for endodontically treating a tooth by orthograde way.

3. Complications
In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

4. Warnings/precautions

- The decision to use an endodontics instrument must be relayed to the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
- Respect recommendations for patients identified as having a risk of infectious endocarditis.
- Contains Nickel and Titanium and must not be used on patients with known allergic sensitivity to these metals.
- Respect the good dental practice in particular by using a dental dam and gloves.
- Use in continuous rotation at the recommended speed.
- In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.
- Use according to the recommended procedure (S7).
- Do not use in endodontics retreatment.
- Respect the expiry date.
- Keep the UDI information which is only on the labelling until the last use.
- Year of manufacture : see labelling.
- Do not use if you doubt on the respect of storage conditions.
- In case of doubt concerning the product identification, do not use.
- Check the condition of the instrument before the use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Inform the manufacturer and the national regulatory authority of any serious incident related to the instrument.
- Be careful to respect the recommended maximum number of uses.
- Instrument supplied non-sterile: follow the reprocessing procedure (S8) before use.

5. Clinical claims
Clinical performance in normal conditions of use :
- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy (respect of the original canal axis)
- Geometry adapted to irrigation
- Maintain of the apical foramen in its original position

Safety performance in normal conditions of use :
- Limited risk of breakage
- Limited risk of shaping aberration
- Limited risk of screwing

6. Characteristics
Type 2 continuous rotating instruments under EN ISO 3630-1:2019 standard:
Maximum recommended torque: see table 'Recommended Torque'.
Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (type 1).

Number of uses : maximum of 5 canals recommended if the instrument is not visually damaged.

Medical device class according to Directive 93/42 and MDR 2017/745: Ila

7. Protocol
EWL : Estimated Working Length

General information on the operating procedure

Start with a catheterisation or exploration of the canal to determine the canal trajectory and the EWL. In the case of difficult canals, passage of the HERO 642 6% provides improved access to the coronary half.

Determine how difficult the channel is so that you can choose the right HERO642 sequence:

- Easy canal : straight and curved <5°
- Canal of medium difficulty: curved 10° < x > 25°
- Difficult canal: curved > 25°

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Easy canals / Blue sequence

1) Conicity of 6%
Set up the coronary access routes, locate the root canal openings and carry out the initial exploration using conventional methods.
On HERO 642® No. 030 6%, adjust the rubber ring to match 1/2 or

2/3 of the EWL.
Penetrate with the instrument in rotation. Use a fast, low amplitude back and forth movement with normal milling pressure.

Successively press in to 1/2 or 2/3 of the depth.

2) Conicity of 4%
On HERO 642® No. 030 4%, adjust the rubber ring to EWL minus 2 mm.

Proceed in the same way as above, combining with circumferential parietal support (surfacing).

If the instrument is no longer progressing, perform a recapitulation with a manual control instrument to check the depth achieved.

3) Conicity of 2%
On HERO 642® No. 030 2%, adjust the rubber ring to the EWL. Carry out a rapid pumping movement to the stop without excessive pressure.

Combine with a circumferential parietal support movement.
Dry before sealing.

Canal of medium difficulty / red sequence :

The crown-down method remains the same, but start with 25 diameter instruments, and follow the red sequence. Use of 5 instruments.

Difficult canal / yellow sequence :
The crown-down method remains the same, but start with 20 diameter instruments, and follow the yellow sequence. Use of 6 instruments.

Enlargement of the apical zone
If the morphology of the canal recommends it (necrotic canal or open apex), it is possible to continue with 2% instruments in diameters 035, 040 and 045 in the apical third.

Note: it is necessary to irrigate abundantly and check apical permeability (using a K-file n010 or n008) between each passage of the instrument and to clean the instrument grooves after each introduction into the canal. Check the instrument between each canal

8. Reprocessing instructions

General recommendations	<ul style="list-style-type: none"> For your own safety, please wear personal protective equipment (gloves, glasses and mask). Do not use cleaning agents containing phenol, aldehyde and alkaline composition. Always respect the instruction for use provided by the product manufacturer.
Limitations on reprocessing	<ul style="list-style-type: none"> Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of medical devices is determined by their function and careful handling. Multiple use of cleaning and re-sterilization cycles may lead to increased risk of file separation. The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. The state of the art and national laws require that validated procedures be followed.
Material needed	<ul style="list-style-type: none"> Gloves, masks, blouse, dental dam, as recommended by the cleaning agents manufacturers. Running Tap water or deionized water. Cleaning agent (neodisher® Septo Active and neodisher® MediZym) Small soft brushes Container Ultrasonic tub or washer-disinfector Class B sterilization apparatus according to EN ISO 17665 <p><i>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i></p>

1	Initial treatment	<ul style="list-style-type: none"> Place used products in a container or a wipe with running tap water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min Rinse the products with running tap water at 20-40 °C for 1 min. If there is time before the next step, make sure the device remains moist by packing it in a wet wipe. Do not exceed 1 hour. <p><i>Remarks:</i> Do not use <i>fixating agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning.</i> Follow the instructions and respect the concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects to devices).</p>	
	Preparation before cleaning	<p>If there are visible impurities on a device, it is recommended to manually pre-clean by brushing under running tap water at 20-40°C for at least 1 min with a soft brush until the impurities have been removed.</p> <p><i>Remark: Follow the instructions and respect the concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects to devices).</i></p>	
	Visual inspection	<p>Inspect the used products and discard the damaged ones (broken, untwisted or anomalously bent).</p>	
2	manual	Cleaning	<p>Insert the products into an ultrasound apparatus beaker.</p> <p>Run ultrasound apparatus during 10-15 min with deionized water and 0,5-2,0% neodisher® MediZym.</p> <p><i>Remarks: Follow the instructions, observe the water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution</i></p>
		Rinsing	<p>Rinse the products with tap water at 20-40°C for 1 min.</p>
		Drying	<p>Dry the products with compressed air until products are visibly dried.</p>
	automatic	Cleaning/Rinsing/Drying	<p>Place the instruments in the tray of the on the slide-in trolley of the washer/disinfector. Ensure that the cleaning agent and water can reach all area of the instrument's blade.</p> <p>Pre-cleaning: cold tap water, 4 min</p> <p>Run cleaning cycle with 0,5%-1,0% neodisher® MediZym at 50°C for 10 min.</p> <p>Final rinsing with deionized water for 3 min.</p> <p>Perform drying.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Follow the instructions and concentrations given by the detergent manufacturer. Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer. Use only approved washer-disinfector according to EN ISO 15883, maintained and validated regularly.
		Visual inspection	<p>Inspect the used products.</p> <p>Repeat steps 4, 5 and 6 if the product is not visibly clean or discard any damaged products.</p>
		Packaging	<p>Place the instruments twice in paper-plastic pouches for steam sterilization in compliance with ISO 11607 and EN 868.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Seal the pouches according to the manufacturer's recommendations. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9	Sterilisation	<p>Sterilize the products with steam vapor</p> <p>Apparatus: class B</p> <p>Temperature: 134°C</p> <p>Time: 18min</p> <p>Minimum drying: 20min</p> <p>Check the physico-chemical indicators and cycle parameters.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> When sterilizing multiple instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded. Place the pouches in the steam sterilizer according to the recommendation given by the sterilizer manufacturer. Use only pre-vacuum air removal steam sterilizer that are matching the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full size sterilizer), with saturated steam. Use only validated sterilizer according to EN ISO 17665, maintained and validated regularly 	
	Storage	<p>Store the products in a dry, clean and dust-free environment at the temperature specified on the paper-plastic pouch by the steam sterilization manufacturer.</p> <p><i>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date). In case of damage, a complete rework should be performed.</i></p>	

9. Storage and transport conditions

10. Disposal
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices..

11. Symbols

	NiTi raw material		Root canal preparation
	Quantity		Assortment
	EN ISO 1797 Type 1		Medical device
	Do not use if the packaging is damaged		Legal manufacturer
	Sterilizable in a steam sterilizer (autoclave) at the specified temperature		Unique Device Identifier

Year of CE-marking : 1998
Revision date : 2024-02-06
Reference : 80000230-A

Micro-Mega SA
12, Rue du Tunnel
25000 Besançon
France
T : +33 3 81 54 42 42
F : +33 3 81 54 42 30
commercial@micro-mega.com
www.micro-mega.com