MicroMega RECI Glider

Reference: 30006063-B

	∅	%	L	Stop	Ring	RPM	Angles	Maximum Recommended Torque (N.cm)
MicroMega RECI Glider	16	03	21 - 25 - 31	Red	White	400 rpm	CCW 170° CW 60°	4.0

Instruction for use	EN
---------------------	----

1. **Indication**
Glide path of root canals during non-surgical endodontic treatment. For use by a dental professional only.

2. **Contraindications**
Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there are no contraindications to the use of MicroMega RECI Glider for endodontically treating a tooth by orthograde route.

3. **Complications**
In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

4. **Warnings/precautions**
• The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
• Patients identified as having a risk of infectious endocarditis.
• Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.
• Respect the good dental practice in particular the use of a dental dam and gloves.

• Use in reciprocating rotation at the recommended angles of rotation.
• In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.
• Use according to the recommended procedure (§7).
• Instruments supplied sterile: Check the integrity of the packaging before use. If damaged, do not use the instrument.
• If applicable, respect the expiry date.
• Keep the UDI information which is only on the labelling until the last use.

• Year of manufacture: see labelling.
• Do not use if doubt on the respect of storage conditions.
• Single use instrument. Do not reuse the instrument. If you do, there is a risk of decreased performance and security characteristics (fracture of the instrument in the tooth).
• In case of doubt concerning the product identification, do not use.
• Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
• Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

Packaging symbol		Sterilized using Ethylene Oxide. Do not re-use
------------------	--	--

5. **Clinical claims**
MicroMega RECI Glider in normal conditions.

Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

Usability:

- Secure progression and limit the interference

Safety:

- Limited risk of breakage
- Limited risk of shaping aberration
- Limited risk of screwing

6. **Characteristics**
Instrument in reciprocating motion in type 4 under EN ISO 3630- 1: 2019 standard: rotation angles are explained in the column «CCW» (Counter-clockwise side) and column «CW» (Clockwise side).

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (Type 1).
Number of uses: single use. 4 canals maximum recommended on the same patient.
Medical device class according to Directive 93/42/EEC and MDR 2017/745: IIa.

7. **Protocol**

General information on the protocol for use:

After opening the access cavity, sodium hypochlorite must be constantly used and renewed before and after each instrument. The sodium hypochlorite must be placed using a syringe in order to be injected as close as possible to the apical region.

Operating procedure :

MicroMega RECI Glider is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.

8. **Reprocessing instructions**

Not applicable.

9. **Storage and transport conditions**

10. **Disposal**
After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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