MicroMega Dual·Pex

INSTRUCTIONS

EN

INSTRUCTIONS

FR

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DE

INSTRUCCIONES

ES

INSTRUCÇÕES

PT

INSTRUCTIES

NL

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RO

KÄYTTÖOHJE

FI

ISTRUZIONI

IT

UTASÍTÁSOK

HU

NÁVOD NA POUŽITIE

SK

BRUKSANVISNING

SV

POKYNY

CS

TALİMATLAR

TR

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EL

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JA

TABLE OF CONTENTS



1. Accessories list (Fig. 1)	3
2. Symbols	3, 132, 133
3. Indication.....	3
4. Contraindications	3
5. Operating and storage conditions.....	3
6. Charging	3
7. Installation	3
8. Function checking procedure	4
9. Volume control	4
10. File location indication	4
11. Use in combination with Dual Move motor.....	5
12. Testing the device.....	5
13. Maintenance.....	5
13. 1. Preamble.....	5
13. 2. Procedure.....	5
14. Main Technical Data	6
15. Electromagnetic Compatibility (EMC) – Declaration of conformity.....	6
16. Warranty.....	8
17. Statement	8

Please carefully read this manual before use, especially warnings.

1. Accessories list (Fig. 1)

1.A	Apex Locator (main unit)	1.B	File clip	1.C	Lip hook	1.D	Measuring wire
1.E	Tester	1.F	Adapter (charger)				

2. Symbols

	Caution		Manufacturer		CE marking		Consult instruction for use
	Product information		Date of manufacture		Direct current		Can be autoclaved up to 135°C
	Serial number		Class II equipment		Comply with WEEE directive		For indoor use only
	Catalogue number		Type BF applied part		Keep away from rain		

3. Indication

This apex locator is used to detect the apex of root canal. This device must only be used in hospital environments, clinics or dental offices by qualified dental employees and not used in the oxygen-rich environment.

4. Contraindications

We recommend not to use the medical device on patients with heart pacemakers or other electrical implants. The device shall not be used to treat children or pregnant women.

5. Operating and storage conditions

- Store and operate the device in a safe and dry environment.
- Avoid the exposure to liquids and direct or indirect heat sources.
- To comply with electromagnetic compatibility (EMC) information, use the device away from fluorescent lamps, radio transmitters, remote controls or mobile RF communication devices.
- If malfunctions occur during the treatment, turn the unit off and contact your after-sales service.
- Never open or repair the device by yourself, otherwise the warranty will not apply.

6. Charging

Fig. 6.A	Connect Dual Pex to the charger and the charger to mains.	Fig. 6.B	Charging is in progress.
-----------------	---	-----------------	--------------------------

- Keep the device away from any heating source and flammable gas.
- Charge the device when low or no power. Too many charge / discharge cycles may reduce the battery shelf-life.
- Use the specified adapter and battery only. Other models may damage the device.
- Do not use the device while charging it.

7. Installation

Fig. 7.A	Insert the electric plug of the measuring wire into the socket on the right-hand side panel of the unit.
-----------------	--

Fig. 7.B	Connect the file clip and the lip hook to the measuring wire.
-----------------	---

- Please plug the measuring wire into the socket in the right direction without forcing. An incorrect connection could compromise the operation and the device accuracy.

8. Function checking procedure

Fig. 8.A	Turn on Dual Pex by touching the power-key. Auto-shut off when unused for 10 minutes.
-----------------	---

Fig. 8.B	Create a short circuit by touching the lip hook with the metal part of the file clip.
-----------------	---

Fig. 8.C	All the horizontal strips shall be displayed on the screen along the full scale and be lighted, a discontinuous beep shall be generated and the "APEX" label shall flash.
-----------------	---

- If the system doesn't meet these indications, check that all the parts are correctly plugged and not damaged.

9. Volume control

Fig. 9	Press the Volume-key several times to adjust the sound volume (alarms and beep when touching keys).
---------------	---

10. File location indication

Connect the lip hook to the patient's lower lip and the file clip to the file that is inserted into the treated root canal. The progression of the file into the canal in reference to the apex will be indicated on the display.

Fig. 10.A	The scale on the display represents the root canal. The green part corresponds to the apical foramen, between the minor apical foramen or apical constriction (marker "1") and the major apical foramen. The red part is over the apex.
------------------	--

Fig. 10.B	The progression of the file towards the apex is showed by horizontal strips displayed along the scale. When the file enters into the front region of the apical foramen, strips are white. At the left of the green scale, the figure indicates the remaining number of strips before reaching the major apical foramen.
------------------	--

Fig. 10.C	Strips become green into the apical zone. The minor and major apical foramens are separated by 7 strips. "00" is displayed when the major apical foramen is reached.
------------------	--

Fig. 10.D	Strips become red over the apex, a rapid beep is emitted and the word "APEX" blinks.
------------------	--

- Open apices, drained canals, oral fluids seepage into the access cavity, root fractures / perforations or canals filled with Gutta Percha may cause erroneous indications.
- Use only the original accessories to guarantee a right functioning and best performances.

- As the green strip number "00" corresponds to the major apical foramen and not the minor one, it is recommended to systematically reduce the working length by 0.5-1 mm.
- Dual Pex is not a measurement device. It indicates that the file is progressively approaching the apex.
- The gingival crevicular fluid, saliva and gingival polyp will interfere with device functioning. Hence it is recommended to correctly isolate the tooth.
- Multi-use accessories (lip hook/ file clip...) which can come in contact with the patient's body during operation must be sterilized by high temperature prior to each use.

11. Use in combination with Dual Move motor

Fig. 11	When Dual Pex is connected to Dual Move, a reference point can be defined to control the file approach. Press the SET-key to define the position of the reference point, between strip 00 and strip 06 ("1" marker). The new position is automatically saved and sent to the connected Dual Move. Please refer to the Dual Move instructions for use for details.
----------------	---

12. Testing the device

It is strongly recommended to regularly control the device using the tester.

Fig. 12	Replace the measuring wire by the tester. The green strips shall appear and the displayed strip number shall be between "02" and "04".
----------------	--

- If the indications are not the expected ones, please check that that the tester is properly connected. If the default remains, stop using Dual Pex and contact your after-sales service.

13. Maintenance

13. 1. Preamble

- Wear gloves and follow carefully the instructions to avoid cross contamination when performing maintenance.
- Refer to your national guidelines, standards and requirements for cleaning, disinfection and sterilization.

- The device and its components have not been disinfected or sterilized before first use. Please follow the necessary steps before the first treatment.
- Cleaning and disinfection protocols to be applied are specific to each component. Be sure to follow the right procedure below when performing daily maintenance.

13. 2. Procedure

- The procedure below shall be followed from step 1 to 5 before the first treatment and after every subsequent treatment.

Operation	Concerned parts	Processing	Warning
1 Cleaning	All	Remove debris under running water (except electronic parts), and then wipe off all components.	
2 Disinfection	Main unit (1.A) Measuring wire (1.D) Mains adapter (1.F)	Wipe the components with a piece of gauze impregnated with Ethanol (Ethanol 70 to 80% vol.). Make sure the components are dry before next step.	If use of washer-disinfectors, use only approved devices according to EN ISO 15883 and refer to their manufacturer's instructions for use.
3 Packing		Pack components in sterilization pouches.	Use pouches which resist to a temperature up to 141°C and comply with EN ISO 11607.

Operation	Concerned parts	Processing	Warning
4 Sterilization	Lip hook (1.C) File clip (1.B)	Steam sterilization at 135°C at least 4 minutes, or at 121°C at least 35 minutes. Minimum drying time after sterilization: 10 minutes.	Use only approved autoclave devices according to EN 13060 or EN 285. The sterilization procedure must comply with ISO 17665. Waiting for cooling before touching.
5 Storage		Keep the components in sterilization packaging in a dry and clean environment.	Check the packaging before using it (packaging integrity, no humidity and validity period), otherwise sterilize again.

- Use only Ethanol for Disinfection (Ethanol 70 to 80% vol.).

14. Main Technical Data

Dimensions - Weight	20 cm x 17.5 cm x 7 cm (package) - 700g
Display	3.5" color LCD
Power supply	Lithium ion battery 3.7V, 1500mAh
Charger	AC 100-240V - 50/60Hz - 0.2VA
Power rating	0.3 W
Degree of protection	IPX 0
Electrical safety class - Applied part	Class II - BF
Operating conditions	Indoor use only 10°C / 40 °C - Relative humidity 30%~75% non-condensing at 0° - Altitude < 2000 m above sea level
Transport and storage conditions	-20 °C / +55 °C - Relative humidity 20% - 80%, non-condensing at > 40 °C - Atmospheric pressure 70 kPa - 106 kPa

15. Electromagnetic Compatibility (EMC) – Declaration of conformity

Guidance and manufacturer's declaration


Dual Pex is intended for use in the electromagnetic environment specified below. The customer or the user of Dual Pex shall make sure it is used in such an environment.

Electromagnetic emission

Emissions test	Compliance	Electromagnetic environment - Guidance
RF emissions - CISPR 11	Group 1	Dual Pex uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions - CISPR 11	Class B	Dual Pex is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/flicker emissions - IEC61000-3-3	Complies	

Electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ±2, 4, 8k, 15V air	±8 kV contact ±2, 4, 8 k 15V air	Floors should be in wood, concrete or ceramic tile. If covered by synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts IEC 61000-4-4	±2 kV 100 kHz repetition frequency	±2 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Line to line: ±0.5kV, ±1 kV	Line to line: ±0.5kV, ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines - IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of devices require continued operation during power mains interruptions, it is recommended that devices be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m 50Hz or 60 Hz	3 A/m 50Hz or 60 Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF - IEC 61000-4-6	3 V 0.15 MHz - 80 MHz, 6V in ISM bands between 0.15 MHz and 80 MHz, 80% AM at 1 KHz	3V 0.15 MHz - 80 MHz, 6V in ISM bands between 0.15 MHz and 80 MHz, 80% AM at 1 KHz	Portable and mobile RF communications equipment should be used no closer to any part of Dual Pex, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 800 MHz ~ 800 MHz $d=2.3\sqrt{P}$ 800 MHz ~ 2.5 GHz
Radiated RF - IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz, 80% AM at 1 KHz	3V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for ratio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicated theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Dual Pex is used exceeds the applicable RF

compliance level above, Dual Pex should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Dual Pex.

** Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the medical device.

Nowadays, many RF wireless equipment's have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The Dual Pex has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and the Dual Pex as recommended below.

Rated maximum output power of transmitter W	Recommended separation distances between portable and mobile RF communications equipment and the Dual Pex.		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.13	0.11	0.22
0.1	0.36	0.37	0.72
1	1.1	1.3	2.1
10	3.5	3.6	7.1
100	13	10	20

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1: at 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

16. Warranty

Inspections and repairs must only be carried out by an approved dealer or by Micro-Mega. Your product is an accurate medical device, produced under maximum quality and test conditions. To ensure an individualized follow-up and traceability of your product, please complete the warranty registration on-line on www.micro-mega.com website within 10 days of your purchase. This will activate the after-sales procedures for your device so that Micro-Mega can give you full satisfaction. Micro-Mega guarantees the battery for one year and the device for two years against any hidden manufacturing defect, excluding faults due to improper use, inappropriate maintenance, accidental damage (fall, impact, etc.), attempted intervention or repairs by personnel not approved by Micro-Mega, attempted product modification.

17. Statement

Non-contractual pictures. Micro-Mega reserves the right to modify the characteristics of its product without prior notice.

COLTENE
MicroMega
Distributed by Micro-Mega
12, rue du Tunnel
25000 Besancon - France

Manufactured by Changzhou Sifary Medical
Technology Co., Ltd.
No. 99, Qingyang Road, Xuejia County, Xinbei
District, Changzhou City, 213000 Jiangsu - China

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Caretechion GmbH