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MicroMega

CE 0459

2Shape / 2Shape Procedure Pack / 2Shape starter pack

Instruction for use

EN

1. Indication

One Flare: flaring during non-surgical endodontics treatment.

2Shape: shaping of root canal during non-surgical endodontics treatment.

One G, K file: glide path of root canals during non-surgical endodontic treatment.

Use by a dental professional only.

2. Contraindications

Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there is no contraindication to the use of 2Shape for endodontically treating a tooth by orthograde.

3. Complications

In cases of complex anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation...) could occur and induce a risk of infectious processes.

4. Warning and precaution

- The decision to use an endodontics instrument must be relayed to the clinical case expertise in particular for a canal anatomy considered as too complex.
- Patients identified as having a risk of infectious endocarditis.
- Contains Nickel and Titanium and should not be used for patient with known allergic sensitivity to these metals.
- Respect the good dental practice in particular the use of dental dam and gloves.
- Use in continuous rotation at the recommended speed.
- In cases of complex anatomy, the maximum torque can be different than the recommended value by the manufacturer.
- Use according to the procedure recommended (§7).
- Do not use in retreatment.
- Instruments supplied non-sterile: Follow the reprocessing instruction (§8) before use.
- Instruments supplied sterile: Check the integrity of the package before use. If damaged, do not use the instrument
- If applicable, respect the expiry date.
- Year of manufacturing: see labelling.
- Do not use if doubt on the respect of storage conditions.
- Instrument single use: Do not reuse the instrument. Otherwise the risk is a decrease of performances and security characteristics (fracture of the instrument in the tooth).
- In case of doubt concerning the product identification, do not use.
- Check the instrument state before the use between each canal. If the instrument is damaged or shows signs of wear, do not use the instrument.
- Inform the manufacturer and the country regulatory authority of any serious incident related to the instrument.

Packaging symbol			Sterile Ethylene Oxide, Do not re-use
Packaging symbol			Not Sterile, Do not re-use
Packaging symbol		-	Not Sterile
Packaging symbol		-	Sterile Ethylene Oxide

5. Clinical claims

2Shape in normal conditions.

Clinical performance:

- Shaping ability: debris removal and cutting efficiency adapted
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintain of the apical foramen in its original position

Usability :

- Improved flexibility through heat treatment

Safety :

- Limited risk of breakage
- Limited risk of aberrations
- No cross-contamination risk

Note: K Files and One G are glide path files used before shaping operation. One Flare is a flaring instrument occurring in the beginning of the shaping.

6. Characteristics

K Files:

- Diameters 10 (purple handle),
- Taper .02 (red stop)
- Material of operative part: Stainless steel
- Lengths: 21 – 25 – 29 mm

Type 2 rotating instruments of EN ISO 3630-1: 2019:

- Maximum recommended torque (not required): see table "Recommended Torque".

	Ø	%	L	Stop	Ring	RPM	Recommended torque (N.cm)
One Flare	25	9	17	-	-	250 400	3
One G	14	3	21 - 25 - 29	●	○	250 400	1.2
2Shape TS1 n25 4%	25	4	21 - 25 - 31	●	●	250 - 400	TS1 : 2N.cm TS2 : 2.5N.cm
2Shape TS2 n25 6%	25	6		●	●		
2Shape F35 n35 4%	35	4	21 - 25 - 31	●	●	250 400	F35 : 2.5N.cm F40 : 2N.cm
2Shape F40 n40 4%	40	4		●	●		


Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (Type 1).

Number of uses: 5 maximum canals recommended.
Medical device class according to 93/42 directive and MDR 2017/745: IIa.

7. Protocol

Protocol 2Shape

Operating procedure :

Progressive movement in three waves (3 up-and-down movements) with upward circumferential filing movement.

Insert the rotating instrument into the root canal until a resistance can be felt. Perform a circumferential brushing movement when feeling the resistance in order to eliminate the primary constraints. Remove the file from the root canal, clean the grooves and irrigate the root canal. Then continue the progressive downward movement. On average, 2 to 3 passages are required before reaching apex.

8. Reprocessing instruction

General recommendations	<ul style="list-style-type: none"> For all metal devices, it is recommended to use anticorrosion disinfecting and cleaning agents. For your own safety, please wear personal protective equipment (gloves, glasses and mask). Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. Always respect the instruction for use provided by the manufacturer of the products.
Limitations on reprocessing	<ul style="list-style-type: none"> Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of the medical devices is determined by their function and careful handling. Multiple use disinfection and re-sterilization cycles may lead to increased risk of file separation. The user must ensure that the processing method used, including re-sources, materials and personnel, is appropriate and meets the applicable requirements. The state of the art and national laws require that validated procedures be followed.
Material needed	<ul style="list-style-type: none"> Gloves, masks, blouse as recommended by the manufacturer of cleaner and detergent Running or deionized water Disinfectant (neodisher® Septo Active) Detergent (neodisher® MediZym) Small soft brushes Container Ultrasonic tub or washer-disinfector Class B sterilization apparatus <p><i>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i></p>

1	<p>Initial Treatment</p> <p>Place used products in a container or a wipe with tap water at 20-40 °C and 1,0 % neodisher® Septo Active during between 5 and 15min.</p> <p>Rinse the products with tap water at 20-40 °C during 1 min.</p> <p>If there is waiting time before next step, make sure the device stays moist by packing it in a wet wipe. Do not exceed 1 hour waiting time.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Do not use fixating agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on devices).
2	<p>Preparation Before Cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C at least 1min with a soft brush until the disappearance of impurities.</p> <p><i>Remark: Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on devices).</i></p>
3	<p>Visual Inspection</p> <p>Inspect the used products and discard the damaged products (broken, untwisted or abnormally bent)</p>

manual	4	<p>Cleaning</p> <p>Insert the products in the ultrasound apparatus beaker</p> <p>Run ultrasound apparatus during 10-30min with tap water and 0,5-2,0% neodisher® MediZym</p> <p><i>Remarks: Follow instructions, observe water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution.</i></p>
	5	<p>Rinsing</p> <p>Rinse the products with tap water at 20-40°C during 1min</p> <p><i>Remark: It is recommended to use deionized water.</i></p>
	6	<p>Drying</p> <p>Dry the products with compressed air until product visibly dried</p>
automatic	4	<p>Cleaning/Rinsing/Drying</p> <p>Place the instruments in the tray of the on the slide-in trolley of the washer/disinfector</p> <p>Perform cleaning cycle with 0,2-1,0% neodisher® MediZym</p> <p>Perform drying</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Disinfection (thermal or chemical-thermal) is not required since the products are sterilized after the cleaning. Follow instructions and concentrations given by the manufacturer of the detergent solution. Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer. The final rinse step should be with deionized water. For other steps follow the water quality defined by the manufacturer. Use only approved washer-disinfector according to EN ISO 15883, maintained and validated regularly.
	7	<p>Visual Inspection</p> <p>Inspect the used products</p> <p>Re-do steps 4-5-6 if the product visibly not clean or discard the damaged products</p>
	8	<p>Packaging</p> <p>Place the instruments in a paper-plastic pouch for steam sterilization compliant to ISO 11607 and EN 868.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent packaging piercing. Seal the pouches according to the recommendation of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9	9	<p>Sterilization</p> <p>Sterilize the products with steam vapor</p> <ul style="list-style-type: none"> Apparatus: class B Minimum temperature: 132°C Minimum time: 3min Absolute pressure: 2,2 bar Minimum drying: 20min <p>Control the physico-chemical indicators and cycle parameters</p> <p><i>Remarks :</i></p> <ul style="list-style-type: none"> The settings temperature at 134 °C and duration of 18 min are recommended for prion inactivation according to French regulation. When sterilizing multiple instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded. Place the pouches in the steam sterilizer according to the recommendation given by the sterilizer manufacturer. Use only pre-vacuum air removal steam sterilizer that are matching the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full size sterilizer), with saturated steam.
	10	<p>Storage</p> <p>Store the products in a dry, clean and dust-free environment at temperature specified by paper-plastic pouch for steam sterilization manufacturer.</p> <p><i>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date). In case of damage, a complete rework should be performed.</i></p>

9. Storage and transport condition

Store sterile instruments in their original packaging, protected from sunlight and heat, in a clean and dry place.

There are no special storage and transport conditions for non-sterile devices.

10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistoury) as per good dentistry practices.

11. Symbols

	"CLASSICS" Handle
	Do not use if packaged is damaged

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Référence : 10300009-A

Micro-Mega SA
12, rue du Tunnel – 25000 BESANCON – FRANCE
Internet : www.micro-mega.com

COLTENE
MicroMega

CE 0459