
MM Bioceramic Putty

CHOOSE YOUR LANGUAGE

ENGLISH.....	5	SLOVENSKY.....	85
ESPAÑOL.....	10	POLSKI.....	91
PORTUGUÊS.....	16	EESTI KEEL.....	96
FRANÇAIS.....	22	LATVIEŠU.....	102
ITALIANO.....	28	LIETUVIŠKAI.....	108
ROMÂNĂ.....	34	ČESKY.....	113
NEDERLANDS.....	40	MAGYAR.....	119
DEUTSCH.....	46	TÜRKÇE.....	125
SVENSKA.....	51	БЪЛГАРСКИ.....	130
NORSK.....	57	СРПСКИ.....	137
DANSKI.....	63	ΕΛΛΗΝΙΚΑ.....	143
SUOMI.....	68	عربي.....	149
HRVATSKI.....	74	MALTI.....	154
SLOVENŠČINA.....	79		

ENGLISH

INTENDED PURPOSE

Ready-to-use bioceramic repair cement, intended for use in patients requiring endodontic procedures, with the function of repairing endodontic complications in the tooth root canal.

INDICATIONS

- A. Treatment of root or furcation perforation via canal.
- B. Treatment of root or furcation perforation via surgery.
- C. Treatment of internal resorption via canal or surgery.
- D. Treatment of external resorption.
- E. Retrograde filling in periapical surgery.
- F. Direct and indirect pulp capping.
- G. Apexification / Apexogenesis
- H. Apexogenesis and Pulpotomy.

INTENDED PATIENT POPULATION

MM Bioceramic Putty should be used in patients with lateral and furcation perforations of the root canal, internal and external resorption, need for retrograde filling, pulp capping, apexification, and apexogenesis.

INTENDED USER

MM Bioceramic Putty is intended to be used by dentists trained to perform endodontic procedures.

INTENDED ENVIRONMENT

Dental offices.

CONTRAINDICATION

- Do not use the product after the expiration date.
- Do not use the product if the packaging is damaged.
- Do not use the product to fill canals, as its flow is inadequate for this purpose.
- Do not use in patients with reported sensitivity to any of the formula's components.

COMPOSITION

Calcium silicate, calcium aluminate, calcium oxide, zirconium oxide, iron oxide, silicon dioxide, and dispersing agent.

TECHNICAL INFORMATION

- Working time: ready-to-use product, no mixing required.

- Setting time: 30 - 120 minutes.

Important: dry with paper points, but do not overdry the cavity. Tissue moisture will initiate the setting reaction of the material.

- High alkalinity.
- High radiopacity.

WARNINGS

- Avoid contact with eyes or skin. In case of contact, rinse immediately with water.
- Avoid contact with oral mucosa. In case of contact, rinse with water and avoid ingestion of the product. If sensitivity persists, seek immediate medical attention.
- Safety and effectiveness of use have not been established in vulnerable populations.
- Keep out of reach of children.

Note: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent health authority where the user and/or patient are established.

PRECAUTIONS

- Use protective glasses, mask, and gloves when handling the product. In case of contact with eyes or skin, wash with water.
- Use rubber dam isolation in all indications for canal use.
- The product is moisture-sensitive; therefore, properly close the package with adequate pressure to prevent drying of the product.
- Before use, check whether the product is hardened at the syringe tip; if so, discard the first portion.
- For endodontic procedures, use the product only after controlling the infection, as material setting does not occur in the acidic pH of infected tissues.
- Cavities where the product will be inserted do not need to be completely dry but must not have excess moisture.

SYRINGE

Before and after use, clean and disinfect the syringe and check whether the silica sachet inside the envelope remains blue. After this procedure, return the syringe properly closed, together with the silica sachet, inside the packaging, which must be correctly closed and stored.

It is recommended to use a disposable plastic cover on the syringe to prevent cross-contamination.

CAPSULE

The MM Bioceramic Putty single-dose capsule is not reusable. After use, discard the capsule and the envelope according to the section "DISPOSAL METHOD."

INSTRUCTIONS FOR USE

CAPSULE

1. Before applying MM Bioceramic Putty single-dose, check if the product shows any damage or defects.
2. Insert the capsule into the gun and remove the capsule cap.
3. Press the gun with slow and uniform pressure to extrude the appropriate amount of material from the capsule.
4. Depending on the application site, MM Bioceramic Putty single-dose may be applied directly with the gun or with an appropriate instrument.
5. Insert the material into the treated area.
6. Complete the treatment according to appropriate methods.
7. Remove the capsule from the gun and discard it.

IMPORTANT:

- Ensure the capsule is at the end of the gun for easier removal. Remove the capsule from the gun and dispose of it after use in accordance with national laws and local authority recommendations for proper disposal.
- Do not reuse the capsule on another patient. Disinfect the gun using the appropriate method before and after use.

A. TREATMENT OF ROOT OR FURCATION PERFORATIONS VIA CANAL

1. Administer anesthesia, place rubber dam isolation, then perform biomechanical preparation of the canal and cleaning of the perforation.
2. Achieve hemostasis and keep the cavity slightly moist.
3. Insert MM Bioceramic Putty and remove excess material in case of overflow.
4. Take radiographs to verify correct placement of the material at the applied site.
5. Fill the remaining canal with gutta-percha cones and a root canal sealer.
6. Perform coronal sealing and restoration.

Important:

- In furcation perforations, apply a layer of glass ionomer as a mechanical barrier, followed by restoration with definitive material.
- Do not place composite resin directly over MM Bioceramic Putty before final set.
- Take radiographs and follow-up for at least two years.

B. TREATMENT OF ROOT OR FURCATION PERFORATIONS VIA SURGERY

1. Administer anesthesia, raise the gingival flap, and perform osteotomy to locate the perforation.

2. Prepare the perforation with an appropriate instrument to facilitate material insertion.
3. Achieve hemostasis and keep the cavity slightly moist.
4. Insert MM Bioceramic Putty into the cavity, adapt with condensers, and remove excess in case of overflow.
5. Reposition the gingival flap and suture.
6. Take radiographs and follow-up for at least two years.

C. TREATMENT OF INTERNAL RESORPTION VIA CANAL OR SURGERY

1. Administer anesthesia and place rubber dam isolation.
2. Remove granulation tissue from the resorption area.
3. Neutralize the environment with calcium hydroxide paste.
4. At the following appointment, remove the calcium hydroxide and perform conventional endodontic treatment up to the resorption level.
5. Apply MM Bioceramic Putty to the resorption site, condensing it against the walls with specific instruments or a slightly moistened sterile cotton pellet.
6. Perform coronal sealing with glass ionomer and restore with material of choice.
7. Take radiographs and follow-up for at least two years.

D. TREATMENT OF EXTERNAL RESORPTION

1. Administer anesthesia, raise the gingival flap, and perform osteotomy to locate the perforation.
2. Remove all affected cementum and dentin, eliminating resorption-inducing cells.
3. Achieve hemostasis, dry the cavity, and condition the affected root surface with citric acid.
4. Insert MM Bioceramic Putty into the slightly moist cavity; adapt the material with condensers and remove excess if extruded.
5. Reposition the gingival flap and suture.
6. Take radiographs and follow-up for at least two years.

Important: In cases of external resorption above the bone crest, MM Bioceramic Putty is not indicated.

E. RETROGRADE FILLING IN PERIAPICAL SURGERY

1. Administer anesthesia, raise the gingival flap, and perform osteotomy to expose the root apex.
2. Resect 3 mm of the root apex containing numerous collateral canals (apical delta).
3. Perform retro-preparation with appropriate instruments.
4. Achieve hemostasis and keep the cavity slightly moist.
5. Insert MM Bioceramic Putty into the retro-cavity, adapt with appropriate condensers, and remove any excess extruded material.

6. Reposition the gingival flap and suture.
7. Take radiographs and follow-up for at least two years.

F. DIRECT AND INDIRECT PULP CAPPING

1. Administer anesthesia and place rubber dam isolation.
 2. Remove caries; if pulp exposure occurs, achieve hemostasis.
 3. Apply MM Bioceramic Putty over the exposed pulp with appropriate instruments, seating it with a slightly moistened sterile cotton pellet.
 4. Cover MM Bioceramic Putty with glass ionomer cement and wait 4 to 6 weeks before placing the final restoration. Glass ionomer can be used as a restoration base.
- Attention: For primary teeth or young permanent teeth with large exposure, perform pulpotomy and follow steps 3 and 4.

G. APEXIFICATION / APEXOGENESIS

First Appointment

1. Administer anesthesia and place rubber dam isolation.
2. Prepare the canal using conventional technique.
3. Apply intracanal medication with calcium hydroxide paste for one week.

Second Appointment

1. Remove the intracanal dressing with files and irrigating solution.
2. Dry the canal with paper points.
3. Create a mechanical barrier by placing a collagen sponge in the apical region to prevent extrusion of MM Bioceramic Putty.
4. Apply MM Bioceramic Putty, condensing it to form a 3–4 mm apical plug.
5. Take radiographs to verify proper canal filling.
6. Place a slightly moistened sterile cotton pellet at the canal entrance followed by temporary restoration.

Third Appointment

1. Remove the temporary restoration and obturate the canal up to the apical plug.
2. Restore with glass ionomer cement.
3. Perform clinical and radiographic follow-up every 3–6 months until radiographic confirmation of hard tissue apical barrier formation.

H. APEXOGENESIS AND PULPOTOMY

1. Administer anesthesia, place rubber dam isolation, then access the pulp chamber.
2. Remove the roof of the pulp chamber and excise coronal pulp using sharp curettes or burs.

3. Section the pulp 0.5 mm below the entrance of the root canal.
4. Thoroughly irrigate the surgical wound with distilled water, saline, or similar solution.
5. Achieve hemostasis and keep the cavity slightly moist.
6. Apply MM Bioceramic Putty over the pulp stump and adapt with a slightly moistened sterile cotton pellet.
7. Take radiographs to verify proper canal filling.
8. Provisionally restore the cavity with glass ionomer.
9. Perform clinical and radiographic follow-up every 3–6 months until root formation is radiographically verified.
10. After root formation is complete, choose between conventional endodontic treatment or only coronal cavity restoration.

STORAGE

- Keep the product in a dry, ventilated place between 15 and 30°C and relative humidity below 60%.
- Do not refrigerate.
- Do not store the product near ammonia, ammonium nitrate, or chlorine-containing products. Avoid using disinfectant solutions containing any of these ingredients.

DISPOSAL METHOD

Dispose of the product in accordance with national laws and recommendations of local authorities for proper disposal. Improper disposal may present biological or physical risks; handle and discard as hazardous biological material if contaminated.

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