

Nickel-Titanium




Ø
1/100 mm

R-Endo®

Ø
20
25
30
35
40

Assortis - Assortied - Assortiert

ISO	ISO	ISO	ISO	ISO	ISO	ISO
Rm	Re	R1	R2	R3	Rs	
4%	12%	8%	6%	4%	4%	




Shank - Manche
Griff

Manual	InGeT®	InGeT®	InGeT®
17 mm	15 mm	15 mm 19 mm 23 mm	25 mm

Packaging - Conditionnement - Packung Packaging - Conditionnement - Packung Packaging - Conditionnement - Packung Packaging - Conditionnement - Packung

Same length, same taper and same number



K Files Limes Feilen		K Reamers Broches Bohrer											Hedstroem Files Racleurs Hedstroemfeilen								
	N°norm.	8	10	15	20	25	30	35	40	45	50	55	60	70	80	90	100	110	120	130	140
	D1 mm	0,08	0,10	0,15	0,20	0,25	0,30	0,35	0,40	0,45	0,50	0,55	0,60	0,70	0,80	0,90	1,00	1,10	1,20	1,30	1,40
	D2 mm	0,14	0,16	0,21	0,26	0,31	0,36	0,41	0,46	0,51	0,56	0,61	0,66	0,76	0,86	0,96	1,06	1,16	1,26	1,36	1,46
	D3 mm	0,40	0,42	0,47	0,52	0,57	0,62	0,67	0,72	0,77	0,82	0,87	0,92	1,02	1,12	1,22	1,32	1,42	1,52	1,62	1,72
	S K Files K Reamers	■			▲					▲									▲		
ISO	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●





MICRO-MEGA®
France Distribution
 5-12, rue du Tunnel
 F-25006 Besançon cedex
 Tel: +33 (0)3 81 54 42 36
 Fax: +33 (0)3 81 54 42 39
 E-mail : commercial@micro-mega.com





MICRO-MEGA®
International distribution
 133 rue de Genève
 CH-1226 Thonex
 Tel: +41 (0)22 860 85 30
 Fax: +41 (0)22 860 85 31
 E-mail : mmid@micro-mega.com

MICRO-MEGA® Dentalvertrieb
GmbH & Co. KG
 Siemensstrasse 5
 D-61267 Neu Anspach
 Tel: +49 (0)6081 4439 0
 Fax: +49 (0)6081 4439 25
 E-mail : info.de@micro-mega.com

www.micro-mega.com



ENDODONTICS The Reference



Nickel-Titanium



Ø
1/100 mm

ENDOFLARE®

HERO 642®

HERO Shaper®

HERO Apical®



ISO

ISO

ISO

ISO

Ø

20

25

30

35

40

45

12%

6%

4%

2%

6%








4%

6%

8%



Assortis - Assortiert - Assortiert



Shank - Manche
Griff

Classics or InGeT®

Classics
(plastic or metal)

Classics or InGeT®

Classics or InGeT®
or Manual

15 mm

21 mm 21 mm 21 mm
25 mm 25 mm 25 mm
29 mm 29 mm 29 mm

21 mm 21 mm
25 mm 25 mm
29 mm 29 mm

25 mm





Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Same length, same taper
and same number



NiTi Systems ● Méthodes NiTi ● Nickel-Titan Methode



InGeT®



Classics




Manual

Stainless Steel



Ø
1/100 mm



	K Files Limes Feilen		K Reamers Broches Bohrer		Hedstroem Files Racleurs Hedstroemfeilen		Heli-Files® (Triocut)
	ISO	MMC	ISO		ISO	MME	ISO
6	●	●	●				
8	●	●	●		●	●	
10	●	●	●		●	●	
15	●	●	●		●	●	●
20	●		●		●	ASS 8-15	●
25	●	ASS 8-15	●		●		●
30	●		●		●		●
35	●		●		●		●
40	●		●		●		●
45	●		●		●		●
50	●		●		●		●
55	●		●		●		●
60	●		●		●		●
70	●		●		●		●
80	●		●		●		●
90	●		●		●		●
100	●		●		●		
110	●		●		●		
120	●		●		●		
130	●		●		●		
140	●		●		●		
	15-40 / 45-80 / 90-140		15-40 / 45-80 / 90-140		15-40 / 45-80 / 90-140		15-40 / 45-80
	21 mm 25 mm 29 mm		21 mm 25 mm 29 mm		21 mm 25 mm 29 mm		21 mm 25 mm 29 mm
	Packaging - Conditionnement - Packung		Packaging - Conditionnement - Packung		Packaging - Conditionnement - Packung		Packaging - Conditionnement - Packung

Assorted
Assortis
Assortiert

Hand use ● Utilisation manuelle ● Handgebrauch

Nerve broaches
Tire-nerfs
Nervnadeln

Stainless Steel



∅
1/100 mm

Pastinject®

Spiral fillers
Bourre-pâte
Wurzelfüller

Engine reamers
Broches Forets
Nervkanalbohrer



MM ISO

ISO

ISO

ISO

-
- ● XXXXF
- ● XXXF
- ● XXF

- ● XF
- ● F
-
- ● M
-
- ● C

-
-

- 6
- 8
- 10
- 15
- 20
- 25
- 30
- 35
- 40
- 45
- 50
- 55
- 60
- 70
- 80
- 90
- 100
- 110
- 120
- 130
- 140

MM 15-80 / ISO 25-60

Assorted
Assortis
Assortiert

25-40

25-40

15-40 / 45-80 / 90-140

- 21 mm
- 29 mm

- 21 mm
- 25 mm
- 29 mm

- 21 mm
- 25 mm
- 29 mm





- 21 mm
- 25 mm
- 29 mm

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung



Contra-angles ● Contre-angles
Winkelstücke

































































Contra-angles
Contre-angles
Winkelstücke
Giromatic®



Stainless Steel



Ø
1/100 mm

	Alésoirs Giromatic® Giromatic® cleansers Giromatic® reibahlen	Giro-Files®	Heli-Giro-Files®	Rispi®
	 	 	 	 
	ISO	ISO	ISO	N° ISO
6				
8				
10				0 
15				1 
20				2 
25				3 
30				4 
35				
40				5 
45				
50				6 
55				
60				
70				
80				
90				
100				
110				
120				
130				
140				
	15-80	15-40 / 45-80 / 90-140	15-40 / 45-80	1-6
	21 mm 25 mm 29 mm	21 mm 25 mm 29 mm	21 mm 25 mm 29 mm	19 mm 22 mm
	Packaging - Conditionnement - Packung	Packaging - Conditionnement - Packung	Packaging - Conditionnement - Packung	Packaging - Conditionnement - Packung
				

Assorted
Assortis
Assortiert

Giromatic®



Stainless Steel





































Ø
1/100 mm

Meca-Shaper®
Shaper-Sonic®
by Dr J.M. Laurichesse

Meca-Rispi®
Rispi-Sonic®
by Dr Riitano & Spina

Finder Spreaders

Finger Pluggers

	 	 	 	 
	ISO	N° ISO	ISO	ISO
6				
8				
10		0 		
15		1 		
20		2 		
25		3 		
30		4 		
35				
40		5 		
45				
50		6 		
55				
60				
70				
80				
90				
100				
110				
120				
130				
140				

Assorted
Assortis
Assortiert

15-40

1-6

15-40

15-40

17-21-25-29 mm
21-25-29 mm Shaper-Sonic®

17-21-25-29 mm
19-22 mm Rispi-Sonic®

21-25 mm





21-25 mm

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung

Packaging - Conditionnement - Packung



Sonic Air®
Mecasonic®



Obturation

