



ROOT CANAL SHAPING

MicroMega One RECI





Exploring opposite directions

micromega.coltene.com

 COLTENE

MicroMega One RECI: single file in reciprocating motion

- **Cutting efficiency:** patented asymmetric cross-section
- **Safety:** cyclic fatigue resistance⁽¹⁾
- **Flexibility:** respect of root canal anatomy thanks to the C.Wire heat treatment⁽²⁾
- **Minimally invasive:** preservation of peri-cervical area and dental hard tissues⁽³⁾



- Single file
- Sterile blister
- 5 sizes (20/.04, 25/.04, 25/.06, 35/.04, 45/.04)
- 3 lengths: > 21 mm
> 25 mm
> 31 mm

MicroMega One RECI, 5 pcs/pack (B5)

60029883	B5 MicroMega One RECI n20 4% L25 Sterile
60029884	B5 MicroMega One RECI n25 4% L25 Sterile
60029885	B5 MicroMega One RECI n25 6% L25 Sterile
60029886	B5 MicroMega One RECI n35 4% L25 Sterile
60029887	B5 MicroMega One RECI n45 4% L25 Sterile
60029899	B5 MicroMega One RECI L25 Assorted Sterile



Complete your treatment procedure with One Flare, MicroMega K-Files, One G, Dual Move & CanalPro Jeni, Gutta Percha - Paper Point



For more information about MicroMega One RECI, please visit our website: micro-mega.com/one-reci



Scan the QR Code to download the detailed list of references or find it on channel.micro-mega.com

- (1) Gianluca Plotino, DDS, PhD,* Hany Mohamed Aly Ahmed, BDS, HDD (Endo), PhD,† Nicola Maria Grande, DDS, PhD,‡ Stephen Cohen, MA, DDS, FICD, and Frederic Bukiet, DDS, MSc, PhD. Current Assessment of Reciprocation in Endodontic Preparation: A Comprehensive Review—Part II: Properties and Effectiveness
- (2) PHILLIP L Tomson, Stéphane R Simon - Contemporary cleaning and shaping of the root canal system - Prim Dent J. 2016;5(2): 46-53
- (3) PJ Van Der Vyver, C. Jonker. Reciprocating instruments in Endodontics: a review of the literature, SADJ. 2014 Oct;69(9):404-9.

©COLTENE

Legal manufacturer:

Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
F +33 3 81 54 42 30
commercial@micro-mega.com
micro-mega.com

MicroMega One RECI | Advertising | Year of CE-marking: 2021
Medical device class IIa according to Directive 93/42/EEC
Notified body: LNE/G-MED | Medical device for dental care,
meant for professional dental use only | See product labelling
and instructions for use.

