


MicroMega **REMOVER**


RETREATMENT



The essential part
of your retreatment

MM MICRO
MEGA

MicroMega REMOVER



- One file to remove the root filling material¹
- Respects the root canal anatomy*
- Efficient without any solvent*
- Safe thanks to non-cutting tip

MicroMega REMOVER

5 instruments, pre-sterilized blister

References	L19	L23
B5 MicroMega Remover n30 7% Classics	20952301	20952311





Diameter

30



Taper

7%



Lengths

19mm

23mm



Wire diameter

1.00mm



Maximum torque

2.5 N.cm



Motor Speed


400rpm

to

800rpm




Continuous rotation



* MicroMega data „Validation report of REMOVER: Pr Franck Diemer, ICRE. Evaluation du One Flare en désobturation. Evaluation d'un instrument de désobturation : le REMOVER.“ 2019.10

¹ Certified for gutta-percha

Visit the product page !



micro-mega.com

Follow us on socials!



@micromegasa



Micro-Mega SA



@micromegasa