

MicroMega REMOVER


RETREATMENT



The essential part
of your retreatment

MM MICRO
MEGA

MicroMega REMOVER



- One file to remove the root filling material¹
- Respects the root canal anatomy*
- Efficient without any solvent*
- Safe thanks to non-cutting tip

MicroMega REMOVER

5 instruments, pre-sterilized blister

References	L19	L23
B5 MicroMega Remover n30 7% Classics	20952301	20952311





Diameter

30



Taper

7%



Lengths

19mm

23mm



Wire diameter

1.00mm



Maximum torque

2.5 N.cm



Motor Speed


400rpm

to

800rpm




Continuous rotation



* MicroMega data „Validation report of REMOVE: Pr Franck Diemer, ICRE. Evaluation du One Flare en désobturation. Evaluation d'un instrument de désobturation : le REMOVE.“ 2019.10

¹ Certified for gutta-percha

Visit the product page !



micro-mega.com

Follow us on socials!



@micromega_international



Micro-Mega SA



Micro-Mega International

60301861-D-03.26

Legal manufacturer :
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega REMOVER
Advertising | Medical device class IIa
according to Directive 93/42/EEC | Notified
body : GMED | Medical device for dental care,
meant for professional dental use only. See
product labelling and instructions for use on
the website : micro-mega.com

CE₀₄₅₉

 **COLTENEGROUP**