

DEBRIDMENT / ROOT CANAL SHAPING

MicroMega

HandFiles



K-Files Hedstroem Files K-Reamers




6 HandFiles per pack

micromega.coltene.com

 **COLTENE**

MicroMega K-Files



K-Files may be used for manual root canal treatments, as well as for mechanized ones, as valuable supplement to NiTi files to achieve the scouting of the root canal with K-Files before the shaping step, K-Files will help doctors determine the working length, with or without the use of an electronic apex locator, also help remove the primary mechanical constraints through an enlargement action.





Available lengths:

21 mm

25 mm

29 mm



Diameters:

006

 to

080



Taper:

2 %


6 instruments per pack, Multiple use, ISO MPL PRIME LINE

Lenght	006	008	010	015	020	025	030	035	040	045	050
L21	20102001L	20102004L	20102007L	20102010L	20102013L	20102016L	20102019L	20102022L	20102025L	20102028L	20102031L
L25	20102002L	20102005L	20102008L	20102011L	20102014L	20102017L	20102020L	20102023L	20102026L	20102029L	20102032L
L29	20102003L	20102006L	20102009L	20102012L	20102015L	20102018L	20102021L	20102024L	20102027L	20102030L	20102033L

Lenght	055	060	070	080
L21	20102034L	20102037L	20102040L	20102043L
L25	20102035L	20102038L	20102041L	20102044L
L29	20102036L	20102039L	20102042L	20102045L

Assortment	008/030	015/040	035/060	045/080
L21	-	20102064L	-	20102067L
L25	20102121L	20102065L	20102122L	20102068L
L29	-	20102066L	-	20102069L

MicroMega Hedstroem Files



Their steeper-angled cutting edges are best suited for vertical dentin cutting, highly efficient enlargement of the root canal and debris removal from the apex to the coronal part, especially in case of irregular canals walls. Last indication relates to removal of separated instruments.



Available lengths:

21 mm

25 mm

29 mm



Diameters:

008

to

080



Taper:

2 %



6 instruments per pack, Multiple use, ISO MPL PRIME LINE

Lenght	008	010	015	020	025	030	035	040	045	050	055
L21	20103010L	20103013L	20103016L	20103019L	20103022L	20103025L	20103028L	20103031L	20103034L	20103037L	20103040L
L25	20103011L	20103014L	20103017L	20103020L	20103023L	20103026L	20103029L	20103032L	20103035L	20103038L	20103041L
L29	20103012L	20103015L	20103018L	20103021L	20103024L	20103027L	20103030L	20103033L	20103036L	20103039L	20103042L

Lenght	060	070	080
L21	20103043L	20103046L	20103049L
L25	20103044L	20103047L	20103050L
L29	20103045L	20103048L	20103051L

Assort- ment	008/030	015/040	045/080
L21	-	-	-
L25	20103108L	20103002L	20103005L
L29	-	-	20103005L

MicroMega K-Reamers



K-Reamers were originally used to find the canals entrances. Today they are mainly used at the end of the root canal preparation as a finisher, to eliminate residual mineral and organic tissues or to place the sealer.



Available lengths:

21 mm

25 mm

29 mm



Diameters:

006

to

080



Taper:

2 %

6 instruments per pack, Multiple use, ISO MPL PRIME LINE

Lenght	006	008	010	015	020	025	030	035	040	045	050
L21	20101001L	20101004L	20101007L	20101010L	20101014L	20101017L	20101020L	20101023L	20101026L	20101029L	20101032L
L25	20101002L	20101005L	20101008L	20101011L	20101015L	20101018L	20101021L	20101024L	20101027L	20101030L	20101033L
L29	20101003L	20101006L	20101009L	20101012L	20101016L	20101019L	20101022L	20101025L	20101028L	20101031L	20101034L

Lenght	055	060	070	080
L21	20101035L	20101038L	20101041L	20101044L
L25	20101036L	20101039L	20101042L	20101045L
L29	20101037L	20101040L	20101043L	20101046L

Assort-ment	015/040	045/080	008/030
L21	20101065L	20101068L	-
L25	20101066L	20101069L	20101117L
L29	20101067L	-	-

© COLTENE – www.coltene.com

Imprint / Publisher

Coltene/Whaledent AG
Feldwiesenstrasse 20
CH-9450 Altstätten
+41 71 757 5300
info.ch@coltene.com

Sales Contact

Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF/UK
+44 1444 235486
info.uk@coltene.com

Service Center

service@coltene.com
00 800 11 0000 12 (UIFN)
Mon.-Fri. 8:00 to 17:00 MESZ

COLTENE