


Shaping

Product upgrade

MicroMega

HERO 642



Root canal shaping sequence
Pack of 6 instruments

micromega.coltene.com

 **COLTENE**

The HERO 642 from COLTENE Micro-Mega is a variable-taper Nickel-Titanium instrument for root canal shaping.

- Continuous rotation at speeds from 300 to 600 rpm
- Crown-down method

Best-sellers

HERO642 N2O 6%



HERO642 N2O 4%



HERO642 N2O 2%



HERO642 N25 4%



HERO642 N30 4%



Product benefits

- The large «residual core» increases resistance to fracture
- Progressive pitch reduces screwing tendency
- Presence of an inactive tip

HERO 642 sequence

- Available lengths :
 - > 21 mm
 - > 25 mm
- Diameters from 20 to 30
- Taper 6, 4 ou 2%

References

20136101	P6 HERO 642 n20 L25 2% Classics	20136114	P6 HERO 642 n20 L21 2% Classics
20136102	P6 HERO 642 n25 L25 2% Classics	20136115	P6 HERO 642 n25 L21 2% Classics
20136103	P6 HERO 642 n30 L25 2% Classics	20136116	P6 HERO 642 n30 L21 2% Classics
20136107	P6 HERO 642 n20 L25 4% Classics	20136120	P6 HERO 642 n20 L21 4% Classics
20136108	P6 HERO 642 n25 L25 4% Classics	20136121	P6 HERO 642 n25 L21 4% Classics
20136109	P6 HERO 642 n30 L25 4% Classics	20136122	P6 HERO 642 n30 L21 4% Classics
20136110	P6 HERO 642 n20 L21 6% Classics	20136132	P6 HERO 642 n20 L25 6% Classics
20136111	P6 HERO 642 n25 L21 6% Classics	20136133	P6 HERO 642 n25 L25 6% Classics
20136112	P6 HERO 642 n30 L21 6% Classics	20136134	P6 HERO 642 n30 L25 6% Classics

Legal manufacturer :
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
F +33 3 81 54 42 30
commercial@micro-mega.com

MicroMega HERO 642
Advertising | Year of CE-marking : 1998 | Medical device class
IIa according to Directive 93/42/EEC | Notified body : GMED |
Medical device for dental care, meant for professional dental
use only. See product labelling and instructions for use on the
website :
micro-mega.com

