

MicroMega **GLIDEPATH**

One Glider
RECI Glider



Creating a path through any obstacle

Product benefits

- Stress decrease and path preparation
- T.Wire : Flexibility
- Flat tip : Excellent guidance and preservation of the root canal anatomy
- Evolving cross section : High resistance and improvement of debris removal
- Available in continuous rotation and reciprocating motion

Ø
Diameter

16



Taper

3%



Available lengths

21mm

25mm

31mm

MicroMega One Glider

5 instruments, Sterile Blister



Continuous rotation



Motor Speed

250 - 400 rpm




Maximum torque

1.2 N.cm

References

B5 MicroMega One Glider n16 3% L21 Ste	60036825
B5 MicroMega One Glider n16 3% L25 Ste	60036827
B5 MicroMega One Glider n16 3% L31 Ste	60036829



Reciprocating angles

170°/60°



Motor Speed

400rpm



Maximum torque


4 N.cm

MicroMega RECI Glider

5 instruments, Sterile Blister

References

B5 MicroMega RECI Glider n16 3% L21 Ste	60036826
B5 MicroMega RECI Glider n16 3% L25 Ste	60036828
B5 MicroMega RECI Glider n16 3% L31 Ste	60036830



80000514-B 06.25

Legal manufacturer:

Micro-Mega SA

12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega One Glider

MicroMega RECI Glider

Advertising | Medical device class IIa according to Directive 93/42/EEC | Notified body : GMED | Medical device for dental care, meant for professional dental use only. See product labelling and instructions for use on the website : micro-mega.com



COLTENE GROUP