

MicroMega

MECHANIZED FILES

# GLIDEPATH

One Glider  
RECI Glider



Creating a path through  
any obstacle

**MM** MICRO  
MEGA

# MicroMega GLIDEPATH

- Stress decrease and path preparation
- T.Wire : Flexibility
- Flat tip : Excellent guidance and preservation of the root canal anatomy
- Evolving cross section : High resistance and improvement of debris removal
- Available in continuous rotation and reciprocating motion



Diameter

16



Taper

3%



Wire diameter

0.8mm



Available lengths

21mm


25mm

31mm

# MicroMega One Glider

5 instruments, pre-sterilized blister

References	
B5 MicroMega One Glider n16 3% L21 Ste	60036825
B5 MicroMega One Glider n16 3% L25 Ste	60036827
B5 MicroMega One Glider n16 3% L31 Ste	60036829



Maximum torque

1.2 N.cm



Motor Speed

250-400 rpm



Continuous rotation



Reciprocating angles

170° / 60°



Motor Speed

400 rpm




Maximum torque

4 N.cm

# MicroMega RECI Glider

5 instruments, pre-sterilized blister

References	
B5 MicroMega RECI Glider n16 3% L21 Ste	60036826
B5 MicroMega RECI Glider n16 3% L25 Ste	60036828
B5 MicroMega RECI Glider n16 3% L31 Ste	60036830



# Visit the product pages !

One Glider



[micro-mega.com](https://micro-mega.com)

RECI Glider



[micro-mega.com](https://micro-mega.com)

## Follow us on socials!



[@micromega\\_international](https://www.instagram.com/micromega_international)



[Micro-Mega SA](https://www.linkedin.com/company/micro-mega-sa)



[Micro-Mega International](https://www.facebook.com/micro-mega-international)

80000514 - D - 03.2026

Legal manufacturer :  
**Micro-Mega SA**  
12 rue du Tunnel  
25000 BESANCON / France  
T +33 3 81 54 42 42  
[mm.sales@coltene.com](mailto:mm.sales@coltene.com)

MicroMega One Glider / MicroMega RECI Glider  
Advertising | Medical device class IIa according  
to Directive 93/42/EEC | Notified body : GMED |  
Medical device for dental care,  
meant for professional dental use only.  
See product labelling and instructions for use on  
the website : [micro-mega.com](https://micro-mega.com)

 0459

 **COLTENEGROUP**