


MicroMega

ROOT CANAL SHAPING

# 2Shape

## Classic & Mini




Sequences of 2 shaping files  
in continuous rotation

**MM** MICRO  
MEGA

# MicroMega 2Shape classic

- Improved debris removal
- Reduced stress
- Increased fracture resistance



Lengths

21mm

25mm

31mm



Wire diameter

1.00mm

1.20mm



Diameters

25

35

40



Tapers

4%

6%



Motor speed

250rpm

to

400rpm



Maximum torque

2.5 N.cm



Continuous rotation

## 2Shape

5 instruments,  
pre-sterilized blister

References	L21	L25	L31
TS1 2Shape N25 4%	20953001	20953011	20953021
TS2 2Shape N25 6%	20953002	20953012	20953022
F35 2Shape N35 6%	20953003	20953013	20953023
F40 2Shape N40 4%	20953004	20953014	20953024

# MicroMega

## 2Shape mini

- Preservation of pericervical dentin
- Elimination of minimal tooth structure
- Respect of root canal anatomy



Maximum torque

2.0 N.cm

2.5 N.cm

TS1 - F40

TS2 - F35



Motor speed

250rpm

to

400rpm



Continuous rotation



Wire diameter

1.00mm



Diameters

25

35


40



Tapers

4%

6%



Lengths

21mm

25mm

31mm

### 2Shape mini

5 instruments,  
pre-sterilized blister

References	L21	L25	L31
TS1 2Shape mini N25 4%	20953201	20953211	20953221
TS2 2Shape mini N25 6%	20953202	20953212	20953222
F35 2Shape mini N35 6%	20953203	20953213	20953223
F40 2Shape mini N40 4%	20953204	20953214	20953224

## 2Shape classic

Guttapercha\*

References	60 pcs.
TS2 2Shape GP 25/.06 SB	60023609
F35 2Shape GP 35/.06 SB	60023610
F40 2Shape GP 40/.04 SB	60023611

Paper points\*

References	144 pcs.
TS2 2Shape PP 25/.06 CP	60023613
F35 2Shape PP 35/.06 CP	60023614
F40 2Shape PP 40/.04 CP	60023615

## 2Shape mini

Guttapercha\*

References	60 pcs.
2Shape mini GP 25/.04 SB	60031543
2Shape mini GP 25/.06 SB	60031545
2Shape mini GP 35/.06 SB	60031546
2Shape mini GP 40/.04 SB	60031544

Paper points\*

References	144 pcs.
2Shape mini PP 25/.04 CP	60031735
2Shape mini PP 25/.06 CP	60031738
2Shape mini PP 35/.06 CP	60031739
2Shape mini PP 40/.04 CP	60031737

\*Products sold by Coltène

# Visit the product page !



[micro-mega.com](https://micro-mega.com)

Legal manufacturer :  
**Micro-Mega SA**  
12 rue du Tunnel  
25000 BESANCON / France  
T +33 3 81 54 42 42  
[mm.sales@coltene.com](mailto:mm.sales@coltene.com)

MicroMega 2Shape / MicroMega 2Shape mini  
Advertising | Years of CE-marking : 2017 (2Shape) et 2021 (2Shape mini) |  
Medical device class IIa according to Directive 93/42/EEC | Notified  
body : GMED | Medical device for dental care, meant for professional  
dental use only. See product labelling and instructions for use on  
the website : [micro-mega.com](https://micro-mega.com)

**CE** 0459

**COLTENEGROUP**