

MOTORS

Cordless endodontic motor

MicroMega

MM 4000



Precision you can trust, performance you can feel

micro-mega.com

MM MICRO
MEGA

3 functions in one device : Integrated apex locator, continuous rotation and reciprocating motion.

- Cordless, for a total freedom of movement
- Intuitive and simplified user interface
- Insulated contra-angle with mini head and slim neck, easy and practical access with clear visibility

ByLedge : Assists in bypassing a ledge

ByPass : Helps navigate a blocked canal or fractured instrument

Pre-programmed Micro-Mega files : + 10 customizable programs



Speed

100-3000 rpm



Torque

0.4-5.0 Ncm



Gear ratio

6:1 - 360° rotation



Reciprocating angle

5° to 400° (adjustable)

References

MM 4000	65003080
---------	----------

Spare parts

Motor Handpiece MM 4000	65003092
Charging base MM 4000	65003093
MM 4000 Box	65003094
Universal AC-Adapter MM 4000	65003121
Black Ring MM 4000	65003122
Battery MM 4000	65003123
Spray Nozzle MM 4000	65003124
File holder MM 4000	65003125
Contra angle MM 4000	65003126
Measuring wire MM 4000	65003127
Lip clip MM 4000	65003128
"O"-ring MM 4000	65003129
Disposable sleeve MM 4000	65003130
Touch probe MM 4000	65003131



80000543-A 05.25

Legal Manufacturer

Guilin Woodpecker

Medical Instruments Co. LTD

Information Industrial Park,
National High-tech Zone.
Guilin

MM 4000

Advertising | Medical device class IIa according to Regulation
MDR 2017/745 | Notified body: TÜV Rheinland | Medical device
for dental care, meant for professional dental use only |
See product labelling and instructions for use on website :
micro-mega.com

COLTENEGROUP