

MicroMega MM4000

CORDLESS ENDO MOTOR



Precision you can trust,
performance you can feel

MM MICRO
MEGA

MicroMega MM 4000

- Cordless, for a total freedom of movement
- Intuitive and simplified user interface
- Insulated contra-angle with mini head and slim neck, easy and practical access with clear visibility

Two new modes :

- **Ledge Bypass** Mode : Assists in bypassing a ledge.
- **Block Bypass** Mode : Helps navigate a blocked canal or a fractured instrument.



Motor Speed

100-3000 rpm



Torque

0.4-5.0 Ncm



Gear ratio

6:1 - 360° rotation



Reciprocating Angles

5° to 400° (adjustable)

MM 4000

References	
MM 4000	65003080
Square parts	
Motor Handpiece MM 4000	65003092
Charging base MM 4000	65003093
MM 4000 Box	65003094
Universal AC-Adapter MM 4000	65003121
Black Ring MM 4000	65003122
Battery MM 4000	65003123
Spray Nozzle MM 4000	65003124
File holder MM 4000	65003125
Contra angle MM 4000	65003126
Measuring wire MM 4000	65003127
Lip clip MM 4000 (X2)	65003128
“O”-ring MM 4000 (X4)	65003129
Disposable sleeve MM 4000 (X10)	65003130
Touch probe MM 4000 (X2)	65003131



Scan for IFU



micro-mega.com

Follow us on socials !



@micromegasa



Micro-Mega SA



@micromegasa

Legal manufacturer :

**Guilin Woodpecker
Medical Instruments Co. LTD**
Information Industrial Park,
National High-tech Zone.
Guilin

MM4000

Advertising | Medical device class IIa
according to Regulation MDR 2017/745 | Notified
body: TÜV Rheinland | Medical device for dental
care, meant for professional dental use only. See
product labelling and instructions for use on
the website : micro-mega.com

CE 0459

 **COLTENE GROUP**