

MicroMega

CORDLESS ENDO MOTOR

# MM4000



Precision you can trust,  
performance you can feel

**MM** MICRO  
MEGA

# MicroMega MM 4000

- Cordless, for a total freedom of movement
- Intuitive and simplified user interface
- Insulated contra-angle with mini head and slim neck, easy and practical access with clear visibility

Two new modes :

- **Ledge Bypass** Mode : Assists in bypassing a ledge.
- **Block Bypass** Mode : Helps navigate a blocked canal or a fractured instrument.



Motor Speed

100-3000 rpm



Torque

0.4-5.0 Ncm



Gear ratio

6:1 - 360° rotation



Reciprocating Angles

5° to 400° (adjustable)

# MM 4000

## References

MM 4000	65003080
---------	----------

## Square parts

Motor Handpiece MM 4000	65003092
-------------------------	----------

Charging base MM 4000	65003093
-----------------------	----------

MM 4000 Box	65003094
-------------	----------

Universal AC-Adapter MM 4000	65003121
------------------------------	----------

Black Ring MM 4000	65003122
--------------------	----------

Battery MM 4000	65003123
-----------------	----------

Spray Nozzle MM 4000	65003124
----------------------	----------

File holder MM 4000	65003125
---------------------	----------

Contra angle MM 4000	65003126
----------------------	----------

Measuring wire MM 4000	65003127
------------------------	----------

Lip clip MM 4000 (X2)	65003128
-----------------------	----------


"O"-ring MM 4000 (X4)	65003129
-----------------------	----------

Disposable sleeve MM 4000 (X10)	65003130
---------------------------------	----------

Touch probe MM 4000 (X2)	65003131
--------------------------	----------



# Scan for IFU



micro-mega.com

# Follow us on socials!



@micromega\_international



Micro-Mega SA



Micro-Mega International

Legal manufacturer :

**Guilin Woodpecker  
Medical Instruments Co. LTD**  
Information Industrial Park,  
National High-tech Zone.  
Guilin

MM4000

Advertising | Medical device class IIa  
according to Regulation MDR 2017/745 | Notified  
body: TÜV Rheinland | Medical device for dental  
care, meant for professional dental use only. See  
product labelling and instructions for use on  
the website : micro-mega.com

