

MicroMega

ROOT CANAL SHAPING

MM TAPER

Scan for IFU



micro-mega.com

Follow us on socials!



@micromegasa



COLTENE
MicroMega



@micromegasa



Shape with Confidence,
Perform with Precision



0459

Legal manufacturer:

Micro-Mega SA

12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega MM TAPER

Advertising | Medical device class IIa
according to Directive 93/42/EEC | Notified
body : GMED | Medical device for dental care,
meant for professional dental use only. See
product labelling and instructions for use on
the website : micro-mega.com

 **COLTENE GROUP**

MM MICRO
MEGA

Product benefits

- Continuous rotation sequence 
- 1 opener, 2 shaping files and 5 finishers
- S.Wire heat treatment: flexibility and resistance to cyclic fatigue

8 instruments, Multiple Use, Sterile Blister

SX 19 / .04v



S1 17 / .02v



S2 20 / .04v



F1 20 / .07v



F2 25 / .08v



F3 30 / .09v



F4 40 / .06v



F5 50 / .06v



Maximum torque

1.5 N.cm	3.1 N.cm	5.1 N.cm
SX - S1	S2 - F1	F2 - F3 - F4 - F5



Motor Speed

250 - 350 rpm



Available lenghts

19mm	21mm	25mm	31mm
SX			

References

B6 MM Taper SX 19 /V04 L19	60040593
B6 MM Taper S1 17/V02 L21	60040594
B6 MM Taper S1 17/V02 L25	60040595
B6 MM Taper S1 17/V02 L31	60040596
B6 MM Taper S2 20/V04 L21	60040597
B6 MM Taper S2 20/V04 L25	60040598
B6 MM Taper S2 20/V04 L31	60040599
B6 MM Taper F1 20/V07 L21	60040600
B6 MM Taper F1 20/V07 L25	60040601
B6 MM Taper F1 20/V07 L31	60040602
B6 MM Taper F2 25/V08 L21	60040603
B6 MM Taper F2 25/V08 L25	60040604

References

B6 MM Taper F2 25/V08 L31	60040605
B6 MM Taper F3 30/V09 L21	60040606
B6 MM Taper F3 30/V09 L25	60040607
B6 MM Taper F3 30/V09 L31	60040608
B6 MM Taper F4 40/V06 L21	60040609
B6 MM Taper F4 40/V06 L25	60040610
B6 MM Taper F4 40/V06 L31	60040611
B6 MM Taper F5 50/V06 L21	60040612
B6 MM Taper F5 50/V06 L25	60040613
B6 MM Taper F5 50/V06 L31	60040614
B6 ASST MM Taper L21	60040615
B6 ASST MM Taper L25	60040616
B6 ASST MM Taper L31	60040617