

MicroMega

ROOT CANAL SHAPING


MM TAPER






Shape with Confidence,
Perform with Precision

MM MICRO
MEGA

Product benefits

- Continuous rotation sequence 
- 1 opener, 2 shaping files and 5 finishers
- S.Wire heat treatment: flexibility and resistance to cyclic fatigue

8 instruments, Multiple Use, Sterile Blister

SX 19 / .04v		
S1 17 / .02v		
S2 20 / .04v		
F1 20 / .07v		
F2 25 / .08v		
F3 30 / .09v		
F4 40 / .06v		
F5 50 / .06v		



Maximum torque

5.1 N.cm

SX - S1

1.5 N.cm

S2 - F1

3.1 N.cm

F2 - F3 - F4 - F5



Motor Speed

250 - 350 rpm



Available lengths

19mm

SX

21mm

25mm

31mm

References

B6 MM Taper SX 19 /V04 L19 **60040593**

B6 MM Taper S1 17/V02 L21 **60040594**

B6 MM Taper S1 17/V02 L25 **60040595**

B6 MM Taper S1 17/V02 L31 **60040596**

B6 MM Taper S2 20/V04 L21 **60040597**

B6 MM Taper S2 20/V04 L25 **60040598**

B6 MM Taper S2 20/V04 L31 **60040599**

B6 MM Taper F1 20/V07 L21 **60040600**

B6 MM Taper F1 20/V07 L25 **60040601**

B6 MM Taper F1 20/V07 L31 **60040602**

B6 MM Taper F2 25/V08 L21 **60040603**

B6 MM Taper F2 25/V08 L25 **60040604**

B6 MM Taper F2 25/V08 L31 **60040605**

References

B6 MM Taper F3 30/V09 L21 **60040606**

B6 MM Taper F3 30/V09 L25 **60040607**

B6 MM Taper F3 30/V09 L31 **60040608**

B6 MM Taper F4 40/V06 L21 **60040609**

B6 MM Taper F4 40/V06 L25 **60040610**

B6 MM Taper F4 40/V06 L31 **60040611**

B6 MM Taper F5 50/V06 L21 **60040612**

B6 MM Taper F5 50/V06 L25 **60040613**

B6 MM Taper F5 50/V06 L31 **60040614**

B6 ASST MM Taper L21 **60040615**

B6 ASST MM Taper L25 **60040616**

B6 ASST MM Taper L31 **60040617**

Scan for IFU



micro-mega.com

Follow us on socials!



@micromega_international



Micro-Mega SA



Micro-Mega International

80000555 - E - 03.26

Legal manufacturer:
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega MM TAPER
Advertising | Medical device class IIa
according to Directive 93/42/EEC | Notified
body : GMED | Medical device for dental care,
meant for professional dental use only. See
product labelling and instructions for use on
the website : micro-mega.com

CE₀₄₅₉

 **COLTENEGROUP**