


MicroMega

ROOT CANAL OBTURATION

# MM Bioceramic Sealer and Putty



Bioceramics :  
Elevating Clinical Standards

**MM** MICRO  
MEGA

# MicroMega Bioceramic

Bioceramic cements **adhere strongly** to dentine, provide **excellent sealing** by limiting microleakage, and **promote healing** by stimulating the formation of new hard tissue, such as dentine or bone.

Developed in accordance with the highest scientific standards, these biomaterials offer a unique combination of excellent sealing ability, exceptional biocompatibility, and easy handling and use.

## Benefits

- Easy application
- Bacteria reduction
- Support healing
- Natural regeneration
- Tight seal

Syringe of 2g

+ 20 tips




### Reference

MM Bioceramic Sealer	65003075
----------------------	----------

### Reference

MM Bioceramic Putty	65003076
---------------------	----------




Syringe of 0.5g

# MM Bioceramic Sealer

MM Bioceramic sealer is used for root canal obturation, combined with a Gutta Percha cone, as a thin layer providing excellent flow and tight seal.

## Product features

- Ready-to-use
- Bioactive
- High Alkalinity (pH 12)
- Excellent flowability
- High radiopacity
- Non staining
- Suitable for cold and warm technique




Working time : 1h  
Setting Time : 2h

# MM Bioceramic Putty

MM Bioceramic Putty is a denser material indicated for root repairs (perforations, resorptions, pulp capping) or retrograde fillings in endodontic surgery.

## Product features

- Ready-to-use
- Bioactive
- High Alkalinity (pH 12)
- Good radiopacity



Visit the product page !



micro-mega.com

Follow us on socials!



@micromegasa



Micro-Mega SA



@micromegasa

80000582 - A - 01.26

Legal manufacturer :

**Angelus**

**Indústria de produtos  
odontológicos S/A**

Waldir Landgraf Street, 101  
Lindoia, Londrina - Brazil

Imported and distributed  
by Micro-Mega SA

MicroMega MM Bioceramic Sealer and Putty  
Advertising | CE 2460 | Medical device class IIa  
according to Regulation MDR 2017/745.

Notified body : DNV (det Norske Veritas).

Medical device for dental care, meant for professional  
dental use only, not reimbursed by Social Security.

See product labelling and instructions  
for use on the website : micro-mega.com

**CE**  
**2460**