


MicroMega

ROOT CANAL SHAPING

2Shape

Classic & Mini




Sequences of 2 shaping files
in continuous rotation

MM MICRO
MEGA

MicroMega 2Shape classic

- Improved debris removal
- Reduced stress
- Increased fracture resistance




2Shape

5 instruments,
pre-sterilized blister

References	L21	L25	L31
TS1 2Shape N25 4 %	20953001	20953011	20953021
TS2 2Shape N25 6 %	20953002	20953012	20953022
F35 2Shape N35 6 %	20953003	20953013	20953023
F40 2Shape N40 4 %	20953004	20953014	20953024

MicroMega 2Shape mini

- Preservation of pericervical dentin
- Elimination of minimal tooth structure
- Respect of root canal anatomy



2Shape mini

5 instruments,
pre-sterilized blister

References	L21	L25	L31
TS1 2Shape mini N25 4 %	20953201	20953211	20953221
TS2 2Shape mini N25 6 %	20953202	20953212	20953222
F35 2Shape mini N35 6 %	20953203	20953213	20953223
F40 2Shape mini N40 4 %	20953204	20953214	20953224

2Shape classic

Guttapercha*

References	60 pcs.
TS2 2Shape GP 25/.06 SB	60023609
F35 2Shape GP 35/.06 SB	60023610
F40 2Shape GP 40/.04 SB	60023611

Paper points*

References	144 pcs.
TS2 2Shape PP 25/.06 CP	60023613
F35 2Shape PP 35/.06 CP	60023614
F40 2Shape PP 40/.04 CP	60023615

2Shape mini

Guttapercha*

References	60 pcs.
2Shape mini GP 25/.04 SB	60031543
2Shape mini GP 25/.06 SB	60031545
2Shape mini GP 35/.06 SB	60031546
2Shape mini GP 40/.04 SB	60031544

Paper points*

References	144 pcs.
2Shape mini PP 25/.04 CP	60031735
2Shape mini PP 25/.06 CP	60031738
2Shape mini PP 35/.06 CP	60031739
2Shape mini PP 40/.04 CP	60031737

*Products sold by Coltène

Visit the product page !



micro-mega.com

Legal manufacturer :
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega 2Shape / MicroMega 2Shape mini
Advertising | Years of CE-marking : 2017 (2Shape) et 2021 (2Shape mini) |
Medical device class IIa according to Directive 93/42/EEC | Notified
body : GMED | Medical device for dental care, meant for professional
dental use only. See product labelling and instructions for use on
the website : micro-mega.com

CE 0459

COLTENE GROUP