

MicroMega

ROOT CANAL SHAPING

# One Curve

## Classic & Mini




Single File : simplicity driving  
performance

**MM** MICRO  
MEGA

# MicroMega One Curve classic

- Resistant to cyclic fatigue
- Suitable for all types of root canals
- Optimized safety
- C-Wire heat treatment



## One Curve

5 instruments,  
presterilized blister

References	L21	L25	L31
One Curve N25 4 %	20949431	20949441	20949451
One Curve N25 6 %	20949401	20949402	20949403
One Curve N35 4 %	20949432	20949442	20949452
One Curve N45 4 %	20949433	20949443	20949453
Assorted Sterile	20949435	20949445	20949455

# MicroMega

# One Curve mini

- Preservation of pericervical dentin
- Elimination of minimal
- Respect of root canal anatomy



Maximum torque

2.5 N.cm



Motor speed

250rpm to 400rpm



Continuous rotation



Wire diameter

1.00mm



Diameters

25

35


40



Tapers

4%

6%



Lengths

21mm

25mm

31mm

## One Curve mini

5 instruments,  
presterilized blister

References	L21	L25	L31
One Curve mini N25 4%	60030017	60030021	60030025
One Curve mini N25 6%	60030018	60030022	60030026
One Curve mini N35 4%	60030019	60030023	60030027
One Curve mini N45 4%	60030020	60030024	60030028
Assorted Sterile	60030033	60030034	60030035

# One Curve classic

Guttapercha\*

Paper points\*

References	60 pcs.
MicroMega OneCurve GP 25/.06 SB	60023612

References	144 pcs.
MicroMega OneCurve PP 25/.06 CP	60023616

# One Curve mini

Guttapercha\*


Paper points\*

References	60 pcs.
One Curve mini GP 25/.04 SB	60031547
One Curve mini GP 25/.06 SB	60031550
One Curve mini GP 35/.04 SB	60031548
One Curve mini GP 45/.04 SB	60037029

References	144 pcs.
One Curve mini PP 25/.04 CP	60031741
One Curve mini PP 25/.06 CP	60031745
One Curve mini PP 35/.04 CP	60031743
One Curve mini PP 45/.04 CP	60031744

\*Products sold by Coltene

## Visit the product page !



[micro-mega.com](https://micro-mega.com)

80000594-A-12.25

**Legal manufacturer :**  
**Micro-Mega SA**  
12 rue du Tunnel  
25000 BESANCON / France  
T +33 3 81 54 42 42  
[mm.sales@coltene.com](mailto:mm.sales@coltene.com)

MicroMega One Curve / MicroMega One Curve mini  
Advertising | Years of CE-marking : 2017 (One Curve) and 2021 (One Curve mini) | Medical device class IIa according to Directive 93/42/EEC | Notified body : GMED | Medical device for dental care, meant for professional dental use only.  
See product labelling and instructions for use on the website : [micro-mega.com](https://micro-mega.com)

CE 0459

 **COLTENEGROUP**