

MicroMega

ROOT CANAL SHAPING

One Curve

Classic & Mini




Single File : simplicity driving
performance

MM MICRO
MEGA

MicroMega

One Curve classic

- Resistant to cyclic fatigue
- Suitable for all types of root canals
- Optimized safety
- C-Wire heat treatment



One Curve

5 instruments,
presterilized blister

References	L21	L25	L31
One Curve N25 4 %	20949431	20949441	20949451
One Curve N25 6 %	20949401	20949402	20949403
One Curve N35 4 %	20949432	20949442	20949452
One Curve N45 4 %	20949433	20949443	20949453
Assorted Sterile	20949435	20949445	20949455

MicroMega

One Curve mini

- Preservation of pericervical dentin
- Elimination of minimal
- Respect of root canal anatomy



Maximum torque

2.0 N.cm

TS1 - F40

2.5 N.cm

TS2 - F35



Motor speed

250rpm

to

400rpm



Continuous rotation



Wire diameter

1.00mm



Diameters

25

35


40



Tapers

4%

6%



Lengths

21mm

25mm

31mm

One Curve mini

5 instruments,
presterilized blister

References	L21	L25	L31
One Curve mini N25 4 %	60030017	60030021	60030025
One Curve mini N25 6 %	60030018	60030022	60030026
One Curve mini N35 4 %	60030019	60030023	60030027
One Curve mini N45 4 %	60030020	60030024	60030028
Assorted Sterile	60030033	60030034	60030035

One Curve classic

Guttapercha*

Paper points*

References	60 pcs.
MicroMega OneCurve GP 25/.06 SB	60023612

References	144 pcs.
MicroMega OneCurve PP 25/.06 CP	60023616

One Curve mini

Guttapercha*

Paper points*

References	60 pcs.
One Curve mini GP 25/.04 SB	60031547
One Curve mini GP 25/.06 SB	60031550
One Curve mini GP 35/.04 SB	60031548
One Curve mini GP 45/.04 SB	60037029

References	144 pcs.
One Curve mini PP 25/.04 CP	60031741
One Curve mini PP 25/.06 CP	60031745
One Curve mini PP 35/.04 CP	60031743
One Curve mini PP 45/.04 CP	60031744

*Products sold by Coltène

Visit the product page !



micro-mega.com

80000594-A-12.25

Legal manufacturer :
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega One Curve / MicroMega One Curve mini
Advertising | Years of CE-marking : 2017 (One Curve) and 2021 (One Curve mini) | Medical device class IIa according to Directive 93/42/EEC | Notified body : GMED | Medical device for dental care, meant for professional dental use only.
See product labelling and instructions for use on the website : micro-mega.com

CE 0459

COLTENEGROUP