


Catalogue HandFiles

MM MICRO
MEGA



2026

CONTENT

MM HandFiles
twisted

5




HandFiles
grinded

7



MM HandFiles
PrimeLine

9




HandFiles


Comparison Overview

	MM K-Files Sterile	MM K-Files	MM K-Files PrimeLine
Manufacturing	Grinded	Twisted	Grinded
Packaging	Sterile	Non Sterile	Non Sterile
Number of Instruments	6	6	6
Sizes	n°6 to n°140	n°6 to n°80	n°6 to n°80
Lengths	21 / 25 / 31 mm	21 / 25 / 31 mm	21 / 25 / 29 mm
Shank	Classic	Classic	PrimeLine
Cutting Efficiency	+	-	+
Flexibility	-	+	-
Resistance to fracture	-	+	-
Also Available	K-Reamers & H-Files	H-Files	K-Reamers & H-Files

Sizes	Flexibility Twisted vs Grinded
n°6	+90%
n°8	+56%
n°10	+89%



Discontinuation
31.12.2026



MicroMega

MM HandFiles

MM K-Files / MM H-Files



Twisted-crafted handfiles for manual scouting or shaping root canals. They serve as an essential complement to Nickel-Titanium instruments, enabling canal location prior to the mechanical shaping stage.

K-Files



↔ Diameter 008 -140

↔ Lengths 21 / 25 / 31 mm

↻ Taper 2%

Pack of 6 instruments



H-Files




↔ Diameter 008 -140


↔ Lengths 21 / 25 / 31 mm

↻ Taper 2%

Pack of 6 instruments



MM K-FILES




Taper

2%

Diameter	L21 mm	L25 mm	L31 mm
006	60040665	60040682	60040701
008	60040666	60040683	60040702
010	60040667	60040684	60040703
015	60040668	60040685	60040704
020	60040669	60040686	60040705
025	60040670	60040687	60040706
030	60040671	60040688	60040707
035	60040672	60040689	60040708
040	60040673	60040690	60040709
045	60040674	60040691	60040710
050	60040675	60040692	60040711
055	60040676	60040693	60040712
060	60040677	60040694	60040713
070	60040678	60040695	60040714
080	60040679	60040696	60040714
Assorted	L21 mm	25 mm	31 mm
008/030	/	60040699	/
015/040	60040680	60040697	60040716
035/060	/	60040700	/
045/080	60040681	60040698	60040717

MM H-FILES



Taper

2%

Diameter	L21 mm	L25 mm	L31 mm
008	60040788	60040804	60040821
010	60040789	60040805	60040822
015	60040790	60040806	60040823
020	60040791	60040807	60040824
025	60040792	60040808	60040825
030	60040793	60040809	60040826
035	60040794	60040810	60040827
040	60040795	60040811	60040828
045	60040796	60040812	60040829
050	60040797	60040813	60040830
055	60040798	60040814	60040831
060	60040799	60040815	60040832
070	60040800	60040816	60040833
080	60040801	60040817	60040834
Assorted	L21 mm	L25 mm	L31 mm
008/030	/	60040818	/
015/040	60040802	60040819	60040835
045/080	60040803	60040820	60040836

MicroMega HandFiles

K-Files / H-Files / K-Reamers



Precision-crafted handfiles for manually scouting or shaping root canals.
Provided in sterile blister.

K-Files



→✕ Diameter 006 -140

↔ Lengths 21 / 25 / 31 mm

↻ Taper 2%

Blister of 6 pre-sterilized instruments



H-Files



→✕ Diameter 008 -140

↔ Lengths 21 / 25 / 31 mm

↻ Taper 2%

Blister of 6 pre-sterilized instruments



K-Reamers



→✕ Diameter 006 -140

↔ Lengths 21 / 25 / 31 mm

↻ Taper 2%

Blister of 6 pre-sterilized instruments




K-Files



Taper

2%


H-Files



Taper

2%

K-Reamers



Taper

2%

Diameter	L21 mm	L25 mm	L31 mm
006	60029717	60029718	60029719
008	60029720	60029721	60029722
010	60029723	60029724	60029725
015	60029726	60029727	60029728
020	60029729	60029730	60029731
025	60029732	60029733	60029734
030	60029735	60029736	60029737
035	60029738	60029739	60029740
040	60029741	60029742	60029743
045	60029744	60029745	60029746
050	60029747	60029748	60029749
055	60029750	60029751	60029752
060	60029753	60029754	60029755
070	60029756	60029757	60029758
080	60029759	60029760	60029761
090	60029762	60029763	60029764
100	60029765	60029766	60029767
110	60029768	60029769	60029770
120	60029771	60029772	60029773
130	60029774	60029775	60029776
140	60029777	60029778	60029779
Assorted	21 mm	25 mm	31 mm
008/030	/	60029992	60029841
015/040	/	/	60029842

Diameter	L21 mm	L25 mm	L31 mm
008	60029780	60029781	60029782
010	60029783	60029784	60029785
015	60029786	60029787	60029788
020	60029789	60029790	60029791
025	60029792	60029793	60029794
030	60029795	60029796	60029797
035	60029798	60029799	60029800
040	60029801	60029802	60029803
045	60029804	60029805	60029806
050	60029807	60029808	60029809
055	60029810	60029811	60029812
060	60029813	60029814	60029815
070	60029816	60029817	60029818
080	60029819	60029820	60029821
090	60029822	60029823	60029824
100	60029825	60029826	60029827
110	60029828	60029829	60029830
120	60029831	60029832	60029833
130	60029834	60029835	60029836
140	60029837	60029838	60029839
Assorted	21 mm	25 mm	31 mm
008/030	/	60029995	60029844
015/040	/	/	60029845

Diameter	L21 mm	L25 mm	L31 mm
006	60029654	60029655	60029656
008	60029657	60029658	60029659
010	60029660	60029661	60029662
015	60029663	60029664	60029665
020	60029666	60029667	60029668
025	60029669	60029670	60029671
030	60029672	60029673	60029674
035	60029675	60029676	60029677
040	60029678	60029679	60029680
045	60029681	60029682	60029683
050	60029684	60029685	60029686
055	60029687	60029688	60029689
060	60029690	60029691	60029692
070	60029693	60029694	60029695
080	60029696	60029697	60029698
090	60029699	60029700	60029701
100	60029702	60029703	60029704
110	60029705	60029706	60029707
120	60029708	60029709	60029710
130	60029711	60029712	60029713
140	60029714	60029715	60029716
Assorted	21 mm	25 mm	31 mm
008/030	/	/	/
015/040	/	60029847	60029848



MicroMega MM HandFiles PrimeLine



MM K-Files / MM H-Files

Switch to MM K-Files / MM H-Files
Available until 12.31.2026

MM K-Files




	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	006	20102001L	20102002L	20102003L
New Ref.		60040665	60040682	60040701
PrimeLine	008	20102004L	20102005L	20102006L
New Ref.		60040666	60040683	60040702
PrimeLine	010	20102007L	20102008L	20102009L
New Ref.		60040667	60040684	60040703
PrimeLine	015	20102010L	20102011L	20102012L
New Ref.		60040668	60040685	60040704
PrimeLine	020	20102013L	20102014L	20102012L
New Ref.		60040669	60040686	60040704
PrimeLine	025	20102016L	20102017L	20102018L
New Ref.		60040670	60040687	60040706
PrimeLine	030	20102019L	20102020L	20102021L
New Ref.		60040671	60040688	60040707
PrimeLine	035	20102022L	20102023L	20102024L
New Ref.		60040672	60040689	60040708
PrimeLine	040	20102025L	20102026L	20102027L
New Ref.		60040673	60040690	60040709
PrimeLine	045	20102028L	20102029L	20102030L
New Ref.		60040674	60040691	60040710

	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	050	20102031L	20102032L	20102033L
New Ref.		60040675	60040692	60040711
PrimeLine	055	20102034L	20102035L	20102036L
New Ref.		60040676	60040693	60040712
PrimeLine	060	20102037L	20102038L	20102039L
New Ref.		60040677	60040694	60040713
PrimeLine	070	20102040L	20102041L	20102042L
New Ref.		60040678	60040695	60040714
PrimeLine	080	20102043L	20102044L	20102045L
New Ref.		60040679	60040696	60040715

Assorted	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	008/030	/	20102121L	/
New Ref.		/	60040699	/
PrimeLine	015/040	20102064L	20102065L	20102066L
New Ref.		60040680	60040697	60040716
PrimeLine	035/060	/	20102122L	/
New Ref.		/	600400700	/
PrimeLine	045/080	20102067L	20102068L	20102069L
New Ref.		60040681	60040698	60040717

MM H-Files



Taper

2%



Lengths

21mm

25mm

31mm




Diameters

008

to

80



	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	008	20103010L	20103011L	20103012L
New Ref.		60040788	60040804	60040821
PrimeLine	010	20103013L	20103014L	20103015L
New Ref.		60040789	60040805	60040822
PrimeLine	015	20103016L	20103017L	20103018L
New Ref.		60040790	60040806	60040823
PrimeLine	020	20103019L	20103020L	20103021L
New Ref.		60040791	60040807	60040834
PrimeLine	025	20103022L	20103023L	20103024L
New Ref.		60040792	60040808	60040825
PrimeLine	030	20103025L	20103026L	20103027L
New Ref.		60040793	60040809	60040826
PrimeLine	035	20103028L	20103029L	20103030L
New Ref.		60040794	60040810	60040827
PrimeLine	040	20103031L	20103032L	20103033L
New Ref.		60040795	60040811	60040828
PrimeLine	045	20103034L	20103035L	20103036L
New Ref.		60040796	60040812	60040829
PrimeLine	050	20103037L	20103038L	20103039L
New Ref.		60040797	60040813	60040830

	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	055	20103040L	20103041L	20103042L
New Ref.		60040798	60040814	60040831
PrimeLine	060	20103043L	20103044L	20103045L
New Ref.		60040799	60040815	60040832
PrimeLine	070	20103046L	20103047L	20103048L
New Ref.		60040800	60040816	60040833
PrimeLine	080	20103049L	20103050L	20103051L
New Ref.		60040801	60040817	60040834

Assorted	Diameter	L21 mm	L25 mm	L29 mm
PrimeLine	008/030	/	20103108L	/
New Ref.		/	60040818	/
PrimeLine	015/040	20103001L	20103002L	20103003L
New Ref.		60040802	60040819	60040835
PrimeLine	035/060	20103004L	20103005L	20103006L
New Ref.		60040803	60040820	60040836
PrimeLine	045/080	20103007L	20103008L	20103009L
New Ref.		/	/	/



Regulatory information

K-Files, K-Reamers, H-Files

CE 0459 | Medical class Isr according to regulation (UE) 2017/745 EEC.

Notified body: GMED.

Legal manufacturer: Micro-Mega SA.

MM K-Files, MM H-Files

Medical device class Iww according to regulation (UE) 2017/745.

Notified body: GMED.


Legal manufacturer: Micro-Mega SA.


PrimeLine MM K-Files, MM H-Files

Medical device class I according to regulation (UE) 2017/745.

Notified body: GMED.

Legal manufacturer: Micro-Mega SA.





Scan to visit our website



micro-mega.com


Follow us on socials



[@micromega_international](https://www.instagram.com/micromega_international)



[Micro-Mega SA](https://www.linkedin.com/company/micro-mega-sa)



[@micromega_international](https://www.facebook.com/micromega_international)

Customer Service Micro-Mega

12 rue du Tunnel
25000 Besançon France
+33 3 81 54 42 36
customer.service@micro-mega.com

80000605-A 05.26

Information products
branded.products.mm@coltene.com

Customer service
mm.sales@coltene.com

Marketing service
mm_marketing@coltene.com

